



High molecular weight gDNA extraction from Western clawed frog (*Xenopus tropicalis*) muscle

19th August 2019

This protocol describes a method to extract high molecular weight genomic DNA from western clawed frog (*Xenopus tropicalis*) muscle, as an example of amphibian tissue. The extraction was performed using the QIAGEN Blood and Cell Culture DNA Midi Kit and part of the genomic DNA was size selected using our protocol "[Size selection of HMW DNA by semi-selective DNA precipitation](#)". Sequencing performance was assessed using the PromethION.

Materials

- 100 mg of frog muscle tissue, stored at -80°C until extraction
- [QIAGEN Blood and Cell Culture DNA Midi Kit](#)
- [QIAGEN TissueRuptor II and probes](#)
- [QIAGEN RNase A](#)
- [Proteinase K](#)
- 2X “size selection buffer” (2.5% w/v PVP 360000 1.2 M NaCl, 20 mM Tris.HCl pH 8)
- [Qubit dsDNA BR Assay Kit \(ThermoFisher Scientific\)](#)
- 70% ethanol in nuclease-free water
- Isopropanol
- TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8)
- 15 ml Falcon tubes
- 1.5 ml Eppendorf DNA LoBind tubes
- Refrigerated centrifuge with capacity for 15 ml Falcon tubes
- Vortex mixer
- Incubator or water bath with capacity for 55°C and agitation capability

Method

● Step 1

Add up to 100 mg of frozen frog muscle tissue to a 15 ml Falcon tube containing 5 ml of buffer G2. Do not allow the tissue to thaw before being placed in the lysis solution.

● Step 2

Homogenise the sample using TissueRuptor II with 2 x 15 second pulses on speed 2. If intact tissue pieces are still visible, repeat until the lysate is homogeneous.

● Step 3

Add 4.5 ml of Buffer G2 and 19 µl of RNase A and mix by inverting the tube.

● Step 4

Incubate the tube at 55 °C for 2 hours, with agitation at 150 rpm.

● Step 5

Equilibrate a QIAGEN Genomic-tip 100/G column with 4 ml of Buffer QBT.

● Step 6

Pour the lysate through the column.

● Step 7

Purify the lysate according to the standard protocol (steps 3–6, pages 50–52).

● Step 8


To maximize DNA yield, we recommend that the elution is performed overnight at room temperature in 200 µl TE buffer.

● Step 9

Take the extracted DNA and perform a size selection using the "[size selection of HMW DNA by semi-selective DNA precipitate](#)" protocol. The expected DNA recovery after size selection is ~25-50%.


Results

- Yield: 50-60 μg
- OD 260/280: 1.98
- OD 260/230: 2.24



Sequencing performance

- Libraries were prepared using the Ligation Sequencing Kit (SQK-LSK109):
- Typical throughput: ★★ (30–60 Gb from FLO-PRO002 flow cells). Throughput from the flow cell may be increased by performing a nuclease wash step at the point where the rate of data acquisition begins to deteriorate due to the accumulation of pores in the “unavailable” or “recovering” state, and then adding a new library.
- Read length profile:



Date	Change note
September 2021	Updated protocol to size select DNA using the size selection of HMW DNA by semi-selective DNA precipitation protocol.

