



High molecular weight gDNA from button mushrooms (*Agaricus bisporus*)

October 2021

This protocol describes a method to extract high molecular weight genomic DNA from button mushrooms (*Agaricus bisporus*) using nuclei isolation, performed as described in [Zhang, M. et al., 2012](#), followed by QIAGEN Genomic-tip purification. Prior to sequencing, 8 µg of genomic DNA was size-selected using the [size selection of HMW DNA by semi-selective DNA precipitation protocol](#). Sequencing performance was determined using the MinION, using the Ligation Sequencing Kit to generate sequencing libraries.

Materials

For nuclei isolation:

- 8 g mushroom gills
- [TissueRuptor II and disposable probes](#)
- [Trizma base](#)
- [KCl](#)
- [Na₂EDTA](#)
- [Spermidine trihydrochloride](#)
- [Spermine tetrahydrochloride](#)
- [NaOH](#)
- [B-mercaptoethanol](#)
- [Triton X-100](#)
- [Sucrose](#)
- ddH₂O
- [Miracloth](#)
- [Cheesecloth](#)
- [Funnel](#)

- 50 ml Falcon tubes
- Refrigerated centrifuge with capacity for 50 ml Falcon tubes
- Ice bucket with ice
- P20, P100, P200 and P1000 pipettes, tips and wide-bore pipette tips

For gDNA extraction and purification:

- [QIAGEN Blood and Cell Culture DNA Midi Kit](#)
- Isopropanol
- 70% ethanol
- [Proteinase K](#)
- Vortex
- 50 ml Falcon tubes
- Refrigerated centrifuge with capacity for 50 ml tubes
- Incubator or water bath with agitation capability and temperature control for 50°C
- TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0)
- P20, P100, P200 and P1000 pipettes, tips and wide-bore pipette tips

Method

Nuclei isolation:

● Step 1

Prepare 10x Homogenisation Buffer (HB) stock as below. Adjust the pH to 9.0-9.4 with NaOH, and store the solution at 4°C.

Tip: When developing this protocol, we made up 1 litre of HB stock and stored it at 4°C for up to a year.

HB stock:

Reagents	Concentration (mM)
Trizma base	100
KCl	800
Na ₂ EDTA	100
Spermidine trihydrochloride	10
Spermine tetrahydrochloride	10

● Step 2

Prepare approximately 90 ml of 1x Homogenisation Buffer (HB) working solution. Keep the solution at 4°C.

HB working solution:

Reagents	Concentration
HB stock	1x
Sucrose	0.5 M

- Step 3

Prepare two 50 ml tubes per sample and add 20 ml of HB working solution to each.

- Step 4

Add 4 g of mushroom gills to each tube.

- Step 5

Homogenise the mushroom sample using the TissueRuptor II at the lowest speed. This will take between 45 seconds to a minute to homogenise.

- Step 6

Add another 20 ml of HB working solution to each tube and invert a few times.

- Step 7

Transfer the tubes to a fume hood and continue the next steps inside the hood until step 20 due to the β -mercaptoethanol toxicity.

- Step 8

Add 0.15% (v/v) of β -mercaptoethanol and invert the tubes 10 times to mix.

- Step 9

Incubate the tubes in a HulaMixer (or equivalent) at 4°C for 10 minutes. If there is no mixer available at 4°C, keep the tubes on ice and gently invert them five times every minute.

- Step 10

Put a funnel on top of two fresh 50 ml Falcon tubes and add a layer of miracloth and two layers of cheesecloth.

● Step 11

Pass the solutions from step 9 through the funnel prepared in step 10. We recommend squeezing the miracloth to maximise the nuclei recovery.

● Step 12

Centrifuge the tubes at 4,000 x g for 20 minutes at 4°C.

● Step 13

Discard the supernatant.

● Step 14

Prepare approximately 45 ml of HB washing solution as below and cool to 4°C.

HB washing solution:

Reagents	Concentration
HB working solution	1x
β-mercaptoethanol	0.15% (v/v)
Triton X-100	0.2% (v/v)

● Step 15

Add 1 ml HB washing solution to the pellets and gently resuspend the pellets by pipetting up and down with a wide-bore pipette tip.

● Step 16

Add 9 ml HB washing solution to both tubes, and gently invert 10 times to mix.

- Step 17

Centrifuge the tubes at 3,100 x g for 15 minutes at 4°C.

- Step 18

Discard the supernatant.

- Step 19

Repeat steps 15-18.

- Step 20

Add 500 µl of 1x HB working solution and gently resuspend the pellet with a wide-bore pipette tip.

DNA extraction:

- Step 1

To each tube with isolated nuclei, add 10 ml of buffer G2 and 100 µl of proteinase K, before inverting the tube 10 times to mix.

- Step 2

Incubate at 50°C for 1 hour with gentle mixing (100 rpm). The tube contents should appear homogeneous but if solid particles are still visible, centrifuge the tubes for a minute at 2,000 x g and transfer the supernatant to fresh tubes.

- Step 3

Equilibrate a Genomic-tip 100/G with 4 ml of buffer QBT.

- Step 4

Pour one of the tubes with lysate into the Genomic-tip 100/G.

- Step 5

Once the content of the first tube has passed through, pour the second tube into the Genomic-tip 100/G.

- Step 6

Purify the lysate according to the [standard QIAGEN protocol](#) (steps 3-5, pages 50-51).

- Step 7

After adding the isopropanol, invert the tubes 10 times and incubate it overnight at -20°C .

- Step 8

Centrifuge the tube at $4,000 \times g$ for 30 minutes at 4°C .

- Step 9

Discard the supernatant and add 5 ml of ice-cold 70% ethanol.

- Step 10

Invert the tube five times and centrifuge at $4,000 \times g$ for 3 minutes at 4°C .

- Step 11

Discard the supernatant and use a clean tissue to dry the walls of the tube.

- Step 12


Elute the pellet in 150 μl TE buffer.

- Optional step

Take 60 μl of the gDNA and size select your sample using the [size selection](#) protocol. About 60% of gDNA is expected to be recovered.

Results

- **Yield:** 15-25 μg
- **A_{260/280}:** 1.92
- **A_{260/230}:** 2.14



Sequencing performance

The library for nanopore sequencing was prepared using the Ligation Sequencing Kit (SQK-LSK110).

- Typical output: ★★★ (>8 Gb in 48 h on FLO-MIN106D) for the Ligation Sequencing Kit. The flow cell was washed using the Flow Cell Wash Kit (EXP-WSH004) and the library re-loaded after ~20 hours of sequencing to maximise flow cell output.
- Read length profile:

