



High molecular weight gDNA extraction from human saliva

19th August 2019

This protocol describes a method to extract high molecular weight genomic DNA from human saliva. The saliva was pooled from multiple individuals to account for variability among samples, and stored in an Isohelix GeneFiX™ Saliva DNA Collection Device. The extraction was carried out using the QIAGEN Blood and Cell Culture DNA Midi Kit followed by size selection with the [size selection of HMW DNA by semi-selective DNA precipitation](#) protocol. Sequencing performance was assessed using the MinION.

Materials

- 1 ml of human saliva
- [Isohelix GeneFiX™ Saliva DNA Collection Device](#)
- [QIAGEN Blood and Cell Culture DNA Midi Kit](#)
- [QIAGEN RNase A](#)
- [Proteinase K](#)
- 2X “size selection buffer” (2.5% w/v PVP 360000 1.2 M NaCl, 20 mM Tris.HCl pH 8)
- [Qubit dsDNA BR Assay Kit \(ThermoFisher Scientific\)](#)
- 70% ethanol in nuclease-free water
- Isopropanol
- TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8)
- 15 ml Falcon tubes
- 1.5 ml Eppendorf DNA LoBind tubes
- Refrigerated centrifuge with capacity for 15 ml Falcon tubes
- Vortex mixer
- Incubator or water bath with capacity for 37°C and 50°C

Method

● Step 1

Collect 1 ml of saliva and store it at room temperature in an Isohelix GeneFiX™ Saliva DNA Collection Device until ready to extract the DNA.

● Step 2

Transfer 1 ml of saliva to a 15 ml Falcon tube and add 19 µl of RNase A.

● Step 3

Add 5 ml of Buffer G2 and 95 µl of Proteinase K to the tube, and mix by inverting the tube.

● Step 4

Incubate the tube for 15 minutes at 37°C, and then for 1 hour at 50°C.

● Step 5

Equilibrate a QIAGEN Genomic-tip 100/G column with 4 ml of Buffer QBT.

● Step 6

Pour the lysate through the column.

● Step 7

Purify the lysate according to the [standard protocol](#) (steps 3–6, pages 50–52).

● Step 8


To maximize DNA yield, we recommend that the elution is performed overnight at room temperature in 150 µl TE buffer.

● Step 9

Take 3 μg of extracted DNA and perform a size selection using the [size selection of HMW DNA by semi-selective DNA precipitation](#) the expected DNA recovery size selection is ~40-50%.

Results


- Yield: 6-12 μg
- OD 260/280: 2.01
- OD 260/230: 2.61



Sequencing performance

Libraries were prepared using the Ligation Sequencing Kit (SQK-LSK109):

- Typical output: ★★★ (8+ Gb in 48 h on FLO-MIN106D) for the Ligation Sequencing Kit, equivalent to the Lambda DNA supplied with the Control Expansion pack (EXP-CTL001).
- Read length profile:



Date	Change note
September 2021	Updated protocol to size select DNA using the size selection of HMW DNA by semi-selective DNA precipitation protocol.
October 2021	Removed reference to QIAGEN TissueRuptor

	II and probes from the materials list.
--	--