

## Materials


- 3 ml of blood, collected in K2-EDTA
- [QIAGEN Puregene Blood Kit](#)
- 70% ethanol in nuclease-free water
- Isopropanol
- TE buffer (1 mM EDTA, pH 8.0)
- 1.5 ml Eppendorf DNA LoBind tubes
- 15 ml Falcon tubes
- Centrifuge and rotor for 15 ml Falcon tubes
- Incubator or water bath (set at 37°C and 50°C)
- Vortex mixer

## Method

1. Perform cell separation and lysis according to the [QIAGEN Puregene Handbook](#) for 3 ml of blood (pages 19–20, steps 1–7):
  - a. Dispense 9 ml of RBC Lysis Solution into a 15 ml Falcon tube.
  - b. Add 3 ml of whole blood and mix by inverting 10 times.
  - c. Incubate for 5 minutes at room temperature. Invert at least once during the incubation.
  - d. Centrifuge for 2 minutes at 2000 x g to pellet the white blood cells.
  - e. Carefully discard the supernatant by pipetting or pouring, leaving approximately 200 µl of the residual liquid and the white blood cell pellet.
  - f. Vortex the tube vigorously to resuspend the pellet in the residual liquid. Vortexing greatly facilitates cell lysis in the next step.
  - g. Add 3 ml of Cell Lysis Solution and pipette mix to lyse the cells or vortex for 10 seconds.
2. **Critical step:** Incubate the samples at 37°C for 30 minutes. If the sample is not homogenous, gently invert the tubes and extend the incubation for another 30 minutes.
3. Purify the lysate according to the [QIAGEN Puregene Handbook](#) for 3 ml of blood (pages 20–21, steps 8–17):
  - a. Add 15 µl of RNase A Solution and mix by inverting 25 times. Incubate for 15 minutes at 37°C. Then incubate for 3 minutes on ice to quickly cool the sample.
  - b. Add 1 ml of Protein Precipitation Solution and vortexed vigorously for 20 seconds at high speed.
  - c. Centrifuge for 5 minutes at 2000 x g. The precipitated proteins should form a tight brown pellet. If the protein pellet is not tight, incubate on ice for 5 minutes and repeat the centrifugation.
  - d. Pipette 3 ml of isopropanol into a clean 15 ml Falcon tube and add the supernatant from the previous step by pouring carefully. Be sure that the protein pellet is not dislodged during pouring.
  - e. Mix by inverting 50 times until the DNA is visible as threads or a clump.
  - f. Centrifuge for 3 minutes at 2000 x g. The DNA may be visible as a small white pellet.
  - g. Carefully discard the supernatant and drain the tube by inverting on a clean piece of absorbant paper, taking care that the pellet remains in the tube.
  - h. Add 3 ml of 70% ethanol and invert several times to wash the DNA pellet.
  - i. Centrifuge for 1 minute at 2000 x g.
  - j. Carefully discard the supernatant. Drain the tube on a clean piece of absorbant paper, taking care that the pellet remains in the tube. Air dry the pellet for 5-10 minutes. The pellet might be loose and easily dislodged. Avoid over-drying the DNA pellet as the DNA will be difficult to dissolve.
  1. **Critical step:** To maximize the DNA yield, we recommend that the elution is performed for 2 hours at 50°C, using 150 µl TE (1 mM EDTA, pH 8.0), occasionally mixing the tube contents by gentle inversion.
  2. Take 3 µg of eluate and perform [aSPRI size selection](#).

## Results


- **Yield:** 20-40  $\mu\text{g}$
- **OD 260/280:** 1.96
- **OD 260/230:** 2.36



## Sequencing performance

Libraries were prepared using the Ligation Sequencing Kit.

- Read length profile:



## Change log

Version	Change
v2, November 2022	Updated Puregene extraction kit name and link
v1, 27th July 2019	Initial protocol publication