

Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 1 ng PolyA+ RNA <input type="checkbox"/> cDNA-PCR Sequencing Kit (SQK-PCS109)	<input type="checkbox"/> Agencourt AMPure XP beads <input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes <input type="checkbox"/> 0.2 ml thin-walled PCR tubes <input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937) <input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water <input type="checkbox"/> 10 mM dNTP solution (e.g. NEB N0447) <input type="checkbox"/> LongAmp Taq 2X Master Mix (e.g. NEB M0287) <input type="checkbox"/> Maxima H Minus Reverse Transcriptase (200 U/μl) with 5x RT Buffer (ThermoFisher, cat # EP0751) <input type="checkbox"/> RNaseOUT™, 40 U/μl (Life Technologies, 10777019) <input type="checkbox"/> Exonuclease I (NEB, M0293)	<input type="checkbox"/> Hula mixer (gentle rotator mixer) <input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes <input type="checkbox"/> Microfuge <input type="checkbox"/> Vortex mixer <input type="checkbox"/> Thermal cycler <input type="checkbox"/> Ice bucket with ice <input type="checkbox"/> Timer <input type="checkbox"/> Pre-chilled freezer block at -20° C for 200 μl tubes (e.g. Eppendorf 022510509) <input type="checkbox"/> Qubit fluorometer (or equivalent for QC check) <input type="checkbox"/> Multichannel pipette capable of 20 - 200 μl <input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
INSTRUCTIONS		NOTES/OBSERVATIONS
Preparing input RNA		
IMPORTANT		
<input type="checkbox"/> Starting from a total RNA sample		
Prepare the RNA in Nuclease-free water <input type="checkbox"/> Transfer 1 ng PolyA+ RNA (or ~50 ng total RNA) into a 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Adjust the volume to up to 9 μl with Nuclease-free water <input type="checkbox"/> Mix by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <input type="checkbox"/> Record the quality, quantity and size of the input RNA.		

cDNA-PCR Sequencing (SQK-PCS109) - PromethION

Version: PPCS_9088_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>IMPORTANT</p> <p>Criteria for input RNA</p> <ul style="list-style-type: none"> <input type="checkbox"/> Average fragment size: ~2 kb <input type="checkbox"/> Input mass, as measured by Qubit RNA HS assay: 1 ng PolyA+ RNA (or ~50 ng total RNA) <input type="checkbox"/> A 260:280 ratio of ~2.0 <input type="checkbox"/> A 260:230 ratio of 2.0-2.2 <input type="checkbox"/> No detergents or surfactants in the buffer 	
<p>Check your flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Click on the Operation wheel icon on the Desktop to load the GUI. <input type="checkbox"/> Load the flow cells to be run into the docking ports of the PromethION. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Wait 10 minutes after loading the flow cells into the PromethION before initiating any experiments or flow cell checks. 	
<ul style="list-style-type: none"> <input type="checkbox"/> On the home screen of the GUI select the flow cells to be run. <input type="checkbox"/> Select the flow cell type from the flow cell drop down. <input type="checkbox"/> Click "Check flow cells" at the bottom left of the screen <input type="checkbox"/> A screen will appear confirming your selections. Click "Start test" to begin the flow cell check 	
<p>When Platform QC completes the message bar will indicate the number of pores found on each flow cell</p>	
<p>Reverse transcription and strand-switching</p>	
<p>Prepare the following reaction in a 0.2 ml PCR tube:</p> <ul style="list-style-type: none"> <input type="checkbox"/> x µl 1 ng PolyA+ RNA (or ~50 ng total RNA) <input type="checkbox"/> 1 µl VNP <input type="checkbox"/> 1 µl 10 mM dNTPs <input type="checkbox"/> 9-x µl RNase-free water <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate at 65° C for 5 minutes and then snap cool on a pre-chilled freezer block. <p>In a separate tube, mix together the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 4 µl 5x RT Buffer <input type="checkbox"/> 1 µl RNaseOUT <input type="checkbox"/> 1 µl Nuclease-free water <input type="checkbox"/> 2 µl Strand-Switching Primer (SSP) <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Add the strand-switching buffer to the snap-cooled, annealed mRNA, mix by flicking the tube and spin down. 	

cDNA-PCR Sequencing (SQK-PCS109) - PromethION

Version: PPCS_9088_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Incubate at 42° C for 2 minutes.</p> <p><input type="checkbox"/> Add 1 µl of Maxima H Minus Reverse Transcriptase. The total volume is now 20 µl.</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Incubate using the following protocol:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Reverse transcription and strand-switching 90 mins @ 42° C (1 cycle) <input type="checkbox"/> Heat inactivation 5 mins @ 85° C (1 cycle) <input type="checkbox"/> Hold @ 4° C 	
<p>Selecting for full-length transcripts by PCR</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> Each PCR reaction uses 5 µl of reverse-transcribed RNA (out of a 20 µl reaction). Therefore, sufficient material is available to perform four PCR reactions per reverse transcription reaction. Do NOT use all 20 µl of the reverse transcription reaction in a single PCR reaction.</p>	
<p>In the same reaction tube, prepare the following reaction at RT:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 25 µl 2x LongAmp Taq Master Mix <input type="checkbox"/> 1.5 µl cDNA Primer (cPRM) <input type="checkbox"/> 18.5 µl Nuclease-free water <input type="checkbox"/> 5 µl Reverse-transcribed RNA sample <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Amplify using the following cycling conditions:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Initial denaturation 30 secs @ 95 °C (1 cycle) <input type="checkbox"/> Denaturation 15 secs @ 95 °C (10-18* cycles) <input type="checkbox"/> Annealing 15 secs @ 62 °C (10-18* cycles) <input type="checkbox"/> Extension 50 secs per kb @ 65 °C (10-18* cycles) <input type="checkbox"/> Final extension 6 mins @ 65 °C (1 cycle) <input type="checkbox"/> Hold @ 4 °C <p><input type="checkbox"/> Add 1 µl of NEB Exonuclease 1 (20 units) directly to each PCR tube.</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate the reaction at 37° C for 15 min, followed by 80° C for 15 min.</p> <p><input type="checkbox"/> Pool the four PCR reactions (total 204 µl) in a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</p> <p><input type="checkbox"/> Add 160 µl of resuspended AMPure XP beads to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p>	

cDNA-PCR Sequencing (SQK-PCS109) - PromethION

Version: PPCS_9088_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 24 µl of Elution Buffer (EB). <input type="checkbox"/> Incubate at RT for 10 mins. <input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 24 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Analyse 1 µl of the amplified DNA for size, quantity and quality. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Sometimes a high-molecular weight product is visible in the wells of the gel when the PCR products are run, instead of the expected smear. These libraries are typically associated with poor sequencing performance. We have found that repeating the PCR with fewer cycles can remedy this. 	
<p>Make up 100-200 fmol of amplified cDNA, and make the volume up to 23 µl in Elution Buffer (EB).</p> <ul style="list-style-type: none"> <input type="checkbox"/> Please check the Mass to Molarity table in the protocol 	
<p>Adapter addition</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Add 1 µl of Rapid Adapter (RAP) to the amplified cDNA library. <input type="checkbox"/> Mix by pipetting and spin down. <input type="checkbox"/> Incubate the reaction for 5 minutes at RT. <input type="checkbox"/> Spin down briefly. 	
<p>The prepared library is used for loading into the flow cell. Store the library on ice until ready to load.</p>	
<p>Priming and loading the flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of PromethION Flush Buffer (PFB) at RT before placing the tubes on ice as soon as thawing is complete. <input type="checkbox"/> Mix the Sequencing Buffer (SQB) and PromethION Flush Buffer (PFB) tubes by vortexing, spin down and return to ice. <input type="checkbox"/> Spin down the Flush Tether (FLT) tube, mix by pipetting, and return to ice. <input type="checkbox"/> Prepare the flow cell priming mix: add 46 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed PromethION Flush Buffer (PFB), and mix by pipetting up and down. 	

cDNA-PCR Sequencing (SQK-PCS109) - PromethION

Version: PPCS_9088_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Load the flow cell(s) into the docking ports within the PromethION</p> <p>Prime the flow cell using the following steps, taking care to avoid the introduction of air bubbles.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Turn the valve to expose the inlet port <input type="checkbox"/> After opening the valve, a small tract of air will be visible beyond the inlet port. Draw back a small volume to remove any air bubbles (a few µls): <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the inlet port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip <input type="checkbox"/> Using a P1000 pipette load 500 µl of the Priming Mix flush into the inlet port of the flow cell, avoiding the introduction of air bubbles <input type="checkbox"/> Wait five minutes. During this time you can prepare your library for loading, as described in the next steps. <input type="checkbox"/> Repeat the priming step with another 500 µl flush <p><input type="checkbox"/> Thoroughly mix the contents of the SQB and LB tubes by pipetting.</p> <p>Prepare the library for loading in a single 1.5 ml Eppendorf DNA LoBind as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 75 µl SQB <input type="checkbox"/> 51 µl LB <input type="checkbox"/> 24 µl DNA library <p>Load your sample</p> <ul style="list-style-type: none"> <input type="checkbox"/> Load 150 µl of your sample through the inlet port <input type="checkbox"/> Close the valve to seal the inlet port <input type="checkbox"/> If no further flow cells are to be loaded, close the PromethION lid <input type="checkbox"/> Wait 10 minutes after loading the flow cells into the PromethION before initiating any experiments 	
<p>Starting a sequencing run</p>	
<p><input type="checkbox"/> Click on the "Operation" wheel icon on the Desktop to load the GUI.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> Wait 10 minutes after loading the flow cells into the PromethION before initiating any experiments or flow cell checks.</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> On the home screen of the GUI select the flow cells to be run. <input type="checkbox"/> Select the flow cell type from the flow cell drop down. <input type="checkbox"/> Click "New experiment" at the bottom left of the screen <input type="checkbox"/> On the New experiment popup screen select the running parameters for your experiment from the individual tabs 	
<p>IMPORTANT</p> <p><input type="checkbox"/> Please do NOT deselect .fast5 in the Output tab. The Dogfish software only passes data to the basecaller after a certain number of files have been written out. Deselecting .fast5 files will result in no basecall data being produced.</p>	

cDNA-PCR Sequencing (SQK-PCS109) - PromethION

Version: PPCS_9088_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>IMPORTANT</p> <p><input type="checkbox"/> Currently basecalling in MinKNOW (through Dogfish) will ignore the options on the new experiment settings screen</p>	
<p><input type="checkbox"/> Once all options have been chosen, click "Start run".</p>	
<p>The home screen will re-appear with your chosen flow cells displaying a progress bar for progression of the sequencing script</p>	
<p>The format and location of your data will depend on the options chosen in the new experiment settings screen:</p>	
<p>Completing a run</p>	
<p>Close down the GUI and power off the PromethION device.</p>	
<p><input type="checkbox"/> Close the GUI window using the close x.</p> <p><input type="checkbox"/> To turn off the PromethION select "Shut Down" from the power options at the top right of the screen.</p> <p><input type="checkbox"/> Once the screen goes blank, turn off the sequencing module via the switch at the rear of the device.</p> <p><input type="checkbox"/> Turn off both Sequencing and Compute Power Supply Units with their respective switches on the back of the units.</p>	
<p>Handling flow cells post-experiment</p>	
<p>If not run for the full 64 hours, flow cells can be stored for subsequent flow cell checks and sequencing experiments.</p> <p>If the flow cell is ready for returning to Oxford Nanopore, please follow the returns and packing instructions.</p> <p>Flongle flow cells are single-use consumables. Once sequencing with a biological sample has occurred and the experiment has finished, dispose of the flow cell via a biological waste procedure.</p>	