

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)

Version: DCB_9091_v109_revT_14Aug2019
Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

Before start checklist

Materials

- 100 ng PolyA+ RNA, or 70-200 ng already-prepared cDNA
- Direct cDNA Sequencing Kit (SQK-DCS109)
- Flow Cell Priming Kit (EXP-FLP002)
- Native Barcoding Expansion 1-12 (EXP-NBD104) and 13-24 (EXP-NBD114) if multiplexing more than 12 samples
- Adapter Mix II Expansion (EXP-AMII001)

Consumables

- Agencourt AMPure XP beads (Beckman Coulter™ cat # A63881)
- NEBNext End repair / dA-tailing Module (E7546)
- NEB Blunt/TA Ligase Master Mix (M0367)
- NEBNext Quick Ligation Module (NEB, E6056)
- 1.5 ml Eppendorf DNA LoBind tubes
- 0.2 ml thin-walled PCR tubes
- Nuclease-free water (e.g. ThermoFisher, AM9937)
- Freshly prepared 70% ethanol in nuclease-free water
- 10 mM dNTP solution (e.g. NEB N0447)
- LongAmp Taq 2X Master Mix (e.g. NEB M0287)
- Maxima H Minus Reverse Transcriptase (200 U/μl) with 5x RT Buffer (ThermoFisher, cat # EP0751)
- RNaseOUT™, 40 U/μl (Life Technologies, cat # 10777019)
- RiboShredder (Epicentre, cat # RS12500), or RNase Cocktail Enzyme Mix (ThermoFisher, cat # AM2286)

Equipment

- Hula mixer (gentle rotator mixer)
- Magnetic rack, suitable for 1.5 ml Eppendorf tubes
- Microfuge
- Vortex mixer
- Thermal cycler
- Ice bucket with ice
- Timer
- Pre-chilled freezer block at -20° C for 200 μl tubes (e.g. Eppendorf cat # 022510509)
- Pipettes and pipette tips P2, P10, P20, P100, P200, P1000

INSTRUCTIONS

NOTES/OBSERVATIONS

Reverse transcription and strand-switching

IMPORTANT

- If you have already prepared your cDNA, use 70-200 ng cDNA and start from the End-prep step.

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)

Version: DCB_9091_v109_revT_14Aug2019
 Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Prepare the RNA in Nuclease-free water</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 100 ng PolyA+ RNA into a 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Adjust the volume to up to 7.5 µl with Nuclease-free water <input type="checkbox"/> Mix by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <p>Prepare the following reaction in a 0.2 ml PCR tube:</p> <ul style="list-style-type: none"> <input type="checkbox"/> x µl poly A+ RNA, 100 ng <input type="checkbox"/> 2.5 µl VNP <input type="checkbox"/> 1 µl 10 mM dNTPs <input type="checkbox"/> 7.5-x µl RNase-free water <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate at 65°C for 5 minutes and then snap cool on a pre-chilled freezer block for 1 minute.</p> <p>In a separate tube, mix together the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 4 µl 5x RT Buffer <input type="checkbox"/> 1 µl RNaseOUT <input type="checkbox"/> 1 µl Nuclease-free water <input type="checkbox"/> 2 µl Strand-Switching Primer (SSP) <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Add the strand-switching buffer to the snap-cooled, annealed mRNA, mix by flicking the tube and spin down.</p> <p><input type="checkbox"/> Incubate at 42°C for 2 minutes in the thermal cycler.</p> <p><input type="checkbox"/> Add 1 µl of Maxima H Minus Reverse Transcriptase. The total volume is now 20 µl.</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Incubate using the following protocol using a thermal cycler:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Reverse transcription and strand-switching 90 mins @ 42°C (1 cycle) <input type="checkbox"/> Heat inactivation 5 mins @ 85°C (1 cycle) <input type="checkbox"/> Hold @ 4°C 	
<p>RNA degradation and second strand synthesis</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Add 1 µl RiboShredder or RNase Cocktail Enzyme Mix (ThermoFisher, cat # AM2286) to the reverse transcription reaction. <input type="checkbox"/> Incubate the reaction for 10 minutes at 37° C in a thermal cycler. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Add 17 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube. 	

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)

Version: DCB_9091_v109_revT_14Aug2019
 Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 20 µl Nuclease-free water.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on the magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 20 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p>Prepare the following reaction in a 0.2 ml thin-walled PCR tube:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 25 µl 2x LongAmp Taq Master Mix <input type="checkbox"/> 2 µl PR2 Primer, 10 µM <input type="checkbox"/> 20 µl Reverse-transcribed sample from above <input type="checkbox"/> 3 µl Nuclease-free water <p>Incubate using the following protocol:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 94 °C <ul style="list-style-type: none"> - Time: 1 mins - No. of cycles: 1 <input type="checkbox"/> 50 °C <ul style="list-style-type: none"> - Time: 1 mins - No. of cycles: 1 <input type="checkbox"/> 65 °C <ul style="list-style-type: none"> - Time: 15 mins - No. of cycles: 1 <input type="checkbox"/> 4 °C <ul style="list-style-type: none"> - Time: ∞ <p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p>	

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)

Version: DCB_9091_v109_revT_14Aug2019
 Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 21 µl Nuclease-free water. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT. <input type="checkbox"/> Pellet the beads on the magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 21 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Analyse 1 µl of the strand-switched DNA for size, quantity and quality. 	
<p>End-prep</p>	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> If you have prepared your own cDNA instead of performing reverse transcription using the Direct cDNA Sequencing Kit, please start this step with 70-200 ng cDNA in 20 µl Nuclease-free water. 	
<p>Perform end repair and dA-tailing of the cDNA sample as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 20 µl cDNA sample <input type="checkbox"/> 30 µl Nuclease-free water <input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer <input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by pipetting and spin down. <input type="checkbox"/> Using a thermal cycler, incubate at 20°C for 5 minutes and 65°C for 5 minutes. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube. <input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. 	

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)



Version: DCB_9091_v109_revT_14Aug2019
 Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 22.5 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute. <input type="checkbox"/> Remove and retain 22.5 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Take forward 22.5 µl of end-prepped cDNA into barcode ligation.</p>	
<p>Barcode ligation</p>	
<p>Add the reagents in the order given below, mixing by flicking the tube between each sequential addition:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 22.5 µl End-prepped DNA <input type="checkbox"/> 2.5 µl Native Barcode <input type="checkbox"/> 25 µl Blunt/TA Ligase Master Mix <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Add 50 µl of resuspended AMPure XP beads to the reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual 70% ethanol. Briefly allow to dry. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 26 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on the magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 26 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)

Version: DCB_9091_v109_revT_14Aug2019
 Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Pool the barcoded samples at the desired ratio to a final volume of 65 µl in a 1.5 ml Eppendorf DNA LoBind tube. Aim for as high a concentration as possible which does not exceed 200 fmol total. If the total volume is >65 µl, perform a 2.5x AMPure clean up and elute in 65 µl of nuclease free water.</p>	
<p>Adapter ligation</p>	
<p>Adapter Mix II Expansion use</p> <p><input type="checkbox"/> Thaw the Wash Buffer (WSB), Elution Buffer (EB) and NEBNext Quick Ligation Reaction Buffer (5x) at RT and mix by vortexing. Then spin down and place on ice. Check the contents of each tube are clear of any precipitate.</p> <p><input type="checkbox"/> Spin down the T4 Ligase and the Adapter Mix II (AMII), and place on ice.</p> <p>Check the contents of each tube are clear of any precipitate and are thoroughly mixed before setting up the reaction.</p> <p><input type="checkbox"/> Check that there is no precipitate present (DTT in the Blunt/TA Master Mix, if used, can sometimes form a precipitate)</p> <p><input type="checkbox"/> Spin down briefly before accurately pipetting the contents into the reaction</p> <p>Taking the pooled and barcoded DNA, perform adapter ligation as follows, mixing by flicking the tube between each sequential addition.</p> <p><input type="checkbox"/> 65 µl 200 fmol pooled barcoded sample</p> <p><input type="checkbox"/> 5 µl Adapter Mix II (AMII)</p> <p><input type="checkbox"/> 20 µl NEBNext Quick Ligation Reaction Buffer (5X)</p> <p><input type="checkbox"/> 10 µl Quick T4 DNA Ligase</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</p> <p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Add 50 µl of resuspended AMPure XP beads to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Place on a magnetic rack, allow beads to pellet and pipette off supernatant.</p> <p><input type="checkbox"/> Add 140 µl of the Wash Buffer (WSB) to the beads. Close the tube lid, and resuspend the beads by flicking the tube. Return the tube to the magnetic rack, allow beads to pellet and pipette off the supernatant.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 13 µl of Elution Buffer (EB).</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on the magnet until the eluate is clear and colourless.</p>	

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)

Version: DCB_9091_v109_revT_14Aug2019
Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Remove and retain 13 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.	
Quantify 1 µl of eluted sample using a Qubit fluorometer.	
The prepared library is used for loading onto the flow cell. Store the library on ice until ready to load.	
Priming and loading the SpotON flow cell	
<input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT before mixing the reagents by vortexing, and spin down the SQB, FB and FLT at RT. <input type="checkbox"/> Prepare the flow cell priming mix: Add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing. <input type="checkbox"/> Open the MinION Mk1B lid and slide the flow cell under the clip. <input type="checkbox"/> Slide the priming port cover clockwise to open the priming port.	
IMPORTANT <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.	
After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles: <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip Note: Visually check that there is continuous buffer from the priming port across the sensor array. <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes. During this time, prepare the library for loading by following the steps below. <input type="checkbox"/> Thoroughly mix the contents of the Loading Beads (LB) by pipetting.	
IMPORTANT <input type="checkbox"/> The Loading Beads (LB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use.	
In a new tube, prepare the library for loading as follows: <ul style="list-style-type: none"> <input type="checkbox"/> 37.5 µl Sequencing Buffer (SQB) <input type="checkbox"/> 25.5 µl Loading Beads (LB), mixed immediately before use <input type="checkbox"/> 12 µl DNA library 	

Direct cDNA sequencing - native barcoding (SQK-DCS109 with EXP-NBD104 and EXP-NBD114)

Version: DCB_9091_v109_revT_14Aug2019
 Last update: 18/05/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of the prepared library to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port and close the priming port. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output. 	
<p>Place the light shield onto the flow cell, as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Carefully place the leading edge of the light shield against the clip. Note: Do not force the light shield underneath the clip. <input type="checkbox"/> Gently lower the light shield onto the flow cell. The light shield should sit around the SpotON cover, covering the entire top section of the flow cell. 	
<p>Close the device lid and set up a sequencing run on MinkNOW.</p>	
<p>Flow cell reuse and returns</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C. <input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol. 	