

Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 100 ng PolyA+ RNA <hr/> <input type="checkbox"/> Direct cDNA Sequencing Kit (SQK-DCS109) <hr/> <input type="checkbox"/> PromethION Flow Cell Priming Kit (EXP-FLP001.PRO.6)	<input type="checkbox"/> Agencourt AMPure XP beads <hr/> <input type="checkbox"/> NEBNext End repair / dA-tailing Module (E7546) <hr/> <input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (M0367) <hr/> <input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes <hr/> <input type="checkbox"/> 0.2 ml thin-walled PCR tubes <hr/> <input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937) <hr/> <input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water <hr/> <input type="checkbox"/> 10 mM dNTP solution (e.g. NEB N0447) <hr/> <input type="checkbox"/> LongAmp Taq 2X Master Mix (e.g. NEB M0287) <hr/> <input type="checkbox"/> Maxima H Minus Reverse Transcriptase (200 U/μl) with 5x RT Buffer (ThermoFisher, cat # EP0751) <hr/> <input type="checkbox"/> RNaseOUT™, 40 U/μl (Life Technologies, 10777019) <hr/> <input type="checkbox"/> RiboShredder (Epicentre, RS12500), or RNase Cocktail Enzyme Mix (ThermoFisher, AM2286)	<input type="checkbox"/> Hula mixer (gentle rotator mixer) <hr/> <input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes <hr/> <input type="checkbox"/> Microfuge <hr/> <input type="checkbox"/> Vortex mixer <hr/> <input type="checkbox"/> Thermal cycler <hr/> <input type="checkbox"/> Ice bucket with ice <hr/> <input type="checkbox"/> Timer <hr/> <input type="checkbox"/> Pre-chilled freezer block at -20° C for 200 μl tubes (e.g. Eppendorf 022510509) <hr/> <input type="checkbox"/> Multichannel pipette capable of 20 - 200 μl <hr/> <input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
INSTRUCTIONS		NOTES/OBSERVATIONS
Preparing input RNA		
Prepare the RNA in Nuclease-free water <input type="checkbox"/> Transfer 100 ng PolyA+ RNA into a DNA LoBind tube <input type="checkbox"/> Adjust the volume to up to 7.5 μl with Nuclease-free water <input type="checkbox"/> Mix by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <input type="checkbox"/> Record the quality, quantity and size of the input RNA.		

Direct cDNA Sequencing (SQK-DCS109) - PromethION

Version: PDCS_9092_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>IMPORTANT</p> <p>Criteria for input RNA</p> <ul style="list-style-type: none"> <input type="checkbox"/> Average fragment size: ~2 kb <input type="checkbox"/> Input mass, as measured by Qubit RNA HS assay: 100 ng <input type="checkbox"/> A 260:280 ratio of ~2.0 <input type="checkbox"/> A 260:230 ratio of 2.0-2.2 <input type="checkbox"/> No detergents or surfactants in the buffer 	
<p>Check your flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Click on the Operation wheel icon on the Desktop to load the GUI. <input type="checkbox"/> Load the flow cells to be run into the docking ports of the PromethION. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Wait 10 minutes after loading the flow cells into the PromethION before initiating any experiments or flow cell checks. 	
<ul style="list-style-type: none"> <input type="checkbox"/> On the home screen of the GUI select the flow cells to be run. <input type="checkbox"/> Select the flow cell type from the flow cell drop down. <input type="checkbox"/> Click "Check flow cells" at the bottom left of the screen <input type="checkbox"/> A screen will appear confirming your selections. Click "Start test" to begin the flow cell check 	
<p>When Platform QC completes the message bar will indicate the number of pores found on each flow cell</p>	
<p>Reverse transcription and strand-switching</p>	
<p>Prepare the following reaction in a 0.2 ml PCR tube:</p> <ul style="list-style-type: none"> <input type="checkbox"/> x µl poly A+ RNA, 100 ng <input type="checkbox"/> 2.5 µl VNP <input type="checkbox"/> 1 µl 10 mM dNTPs <input type="checkbox"/> 7.5-x µl RNase-free water <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate at 65° C for 5 minutes and then snap cool on a pre-chilled freezer block. <p>In a separate tube, mix together the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 4 µl 5x RT Buffer <input type="checkbox"/> 1 µl RNaseOUT <input type="checkbox"/> 1 µl Nuclease-free water <input type="checkbox"/> 2 µl Strand-Switching Primer (SSP) <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Add the strand-switching buffer to the snap-cooled, annealed mRNA, mix by flicking the tube and spin down. 	

Direct cDNA Sequencing (SQK-DCS109) - PromethION

Version: PDCS_9092_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Incubate at 42° C for 2 minutes.</p> <p><input type="checkbox"/> Add 1 µl of Maxima H Minus Reverse Transcriptase. The total volume is now 20 µl.</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Incubate using the following protocol:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Reverse transcription and strand-switching 90 mins @ 42° C (1 cycle) <input type="checkbox"/> Heat inactivation 5 mins @ 85° C (1 cycle) <input type="checkbox"/> Hold @ 4° C 	
<p>RNA degradation and second strand synthesis</p>	
<p><input type="checkbox"/> Add 1 µl RiboShredder or RNase Cocktail Enzyme Mix (ThermoFisher, AM2286) to the reverse transcription reaction.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at 37° C.</p> <p><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</p> <p><input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube.</p> <p><input type="checkbox"/> Add 17 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 20 µl Nuclease-free water.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p> <p><input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 20 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p>Prepare the following reaction in a 0.2 ml thin-walled PCR tube:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 25 µl 2x LongAmp Taq Master Mix <input type="checkbox"/> 2 µl PR2 Primer (PR2) <input type="checkbox"/> 20 µl Reverse-transcribed sample from above <input type="checkbox"/> 3 µl Nuclease-free water 	

Direct cDNA Sequencing (SQK-DCS109) - PromethION

Version: PDCS_9092_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Incubate using the following protocol:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 94 °C 1 mins 1 <input type="checkbox"/> 50 °C 1 mins 1 <input type="checkbox"/> 65 °C 15 mins 1 <input type="checkbox"/> 4 °C ∞ <ul style="list-style-type: none"> <input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing. <input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube. <input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 21 µl Nuclease-free water. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT. <input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 21 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Analyse 1 µl of the strand-switched DNA for size, quantity and quality. 	
<p>End-prep</p>	
<p>Perform end repair and dA-tailing of fragmented DNA as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 20 µl cDNA sample <input type="checkbox"/> 30 µl Nuclease-free water <input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer <input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by pipetting and spin down. <input type="checkbox"/> Using a thermal cycler, incubate at 20° C for 5 minutes and 65° C for 5 mins. <input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing. <input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube. 	

Direct cDNA Sequencing (SQK-DCS109) - PromethION

Version: PDCS_9092_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 30 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 30 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Take forward 30 µl of end-prepped cDNA into adapter ligation.</p>	
<p>Adapter ligation</p>	
<p>Check the contents of each tube are clear of any precipitate and are thoroughly mixed before setting up the reaction.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Mix the contents of each tube by flicking <input type="checkbox"/> Check that there is no precipitate present (DTT in the Blunt/TA Master Mix can sometimes form a precipitate) <input type="checkbox"/> Spin down briefly before accurately pipetting the contents in the reaction <p>Taking the end-prepped DNA, perform adapter ligation as follows, mixing by flicking the tube between each sequential addition.</p> <ul style="list-style-type: none"> <input type="checkbox"/> 30 µl End-prepped DNA <input type="checkbox"/> 5 µl Adapter Mix <input type="checkbox"/> 50 µl Blunt/TA Ligation Master Mix <input type="checkbox"/> 15 µl Nuclease-free water <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. 	
<p>Adapted and tethered DNA library.</p>	
<p>AMPure XP bead binding</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing. 	

Direct cDNA Sequencing (SQK-DCS109) - PromethION

Version: PDCS_9092_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the adapter ligation reaction from the previous step and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Place on magnetic rack, allow beads to pellet and pipette off supernatant. <input type="checkbox"/> Add 200 µl of ABB Buffer (ABB) to the beads. Close the tube lid, and resuspend the beads by pipetting. Return the tube to the magnetic rack, allow beads to pellet and pipette off the supernatant. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 25 µl of Elution Buffer (EB). <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 25 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Quantify 1 µl of eluted cDNA using a Qubit fluorometer - recovery aim ~60 fmol. 	
<p>IMPORTANT</p> <p>Please be aware that the flow cell's pore occupancy could be compromised when loading lower amounts of cDNA. Please use the table below as a guide:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Please check the Mass to Molarity table in the protocol 	
<p>The prepared library is used for loading into the flow cell. Store the library on ice until ready to load.</p>	
<p>Priming and loading the flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of PromethION Flush Buffer (PFB) at RT before placing the tubes on ice as soon as thawing is complete. <input type="checkbox"/> Mix the Sequencing Buffer (SQB) and PromethION Flush Buffer (PFB) tubes by vortexing, spin down and return to ice. <input type="checkbox"/> Spin down the Flush Tether (FLT) tube, mix by pipetting, and return to ice. <input type="checkbox"/> Prepare the flow cell priming mix: add 46 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed PromethION Flush Buffer (PFB), and mix by pipetting up and down. <input type="checkbox"/> Load the flow cell(s) into the docking ports within the PromethION 	

Direct cDNA Sequencing (SQK-DCS109) - PromethION

Version: PDCS_9092_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Prime the flow cell using the following steps, taking care to avoid the introduction of air bubbles.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Turn the valve to expose the inlet port <input type="checkbox"/> After opening the valve, a small tract of air will be visible beyond the inlet port. Draw back a small volume to remove any air bubbles (a few µls): <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the inlet port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip <input type="checkbox"/> Using a P1000 pipette load 500 µl of the Priming Mix flush into the inlet port of the flow cell, avoiding the introduction of air bubbles <input type="checkbox"/> Wait five minutes. During this time you can prepare your library for loading, as described in the next steps. <input type="checkbox"/> Repeat the priming step with another 500 µl flush <p><input type="checkbox"/> Thoroughly mix the contents of the SQB and LB tubes by pipetting.</p> <p>Prepare the library for loading in a single 1.5 ml Eppendorf DNA LoBind as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 75 µl SQB <input type="checkbox"/> 51 µl LB <input type="checkbox"/> 24 µl DNA library <p>Load your sample</p> <ul style="list-style-type: none"> <input type="checkbox"/> Load 150 µl of your sample through the inlet port <input type="checkbox"/> Close the valve to seal the inlet port <input type="checkbox"/> If no further flow cells are to be loaded, close the PromethION lid <input type="checkbox"/> Wait 10 minutes after loading the flow cells into the PromethION before initiating any experiments 	
<p>Starting a sequencing run</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Click on the "Operation" wheel icon on the Desktop to load the GUI. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Wait 10 minutes after loading the flow cells into the PromethION before initiating any experiments or flow cell checks. 	
<ul style="list-style-type: none"> <input type="checkbox"/> On the home screen of the GUI select the flow cells to be run. <input type="checkbox"/> Select the flow cell type from the flow cell drop down. <input type="checkbox"/> Click "New experiment" at the bottom left of the screen <input type="checkbox"/> On the New experiment popup screen select the running parameters for your experiment from the individual tabs 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Please do NOT deselect .fast5 in the Output tab. The Dogfish software only passes data to the basecaller after a certain number of files have been written out. Deselecting .fast5 files will result in no basecall data being produced. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Currently basecalling in MinKNOW (through Dogfish) will ignore the options on the new experiment settings screen 	

Direct cDNA Sequencing (SQK-DCS109) - PromethION

Version: PDCS_9092_v109_revA_04Feb2019
 Last update: 04/02/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Once all options have been chosen, click "Start run".	
The home screen will re-appear with your chosen flow cells displaying a progress bar for progression of the sequencing script	
The format and location of your data will depend on the options chosen in the new experiment settings screen:	
Completing a run	
Close down the GUI and power off the PromethION device.	
<input type="checkbox"/> Close the GUI window using the close x. <input type="checkbox"/> To turn off the PromethION select "Shut Down" from the power options at the top right of the screen. <input type="checkbox"/> Once the screen goes blank, turn off the sequencing module via the switch at the rear of the device. <input type="checkbox"/> Turn off both Sequencing and Compute Power Supply Units with their respective switches on the back of the units.	
Handling flow cells post-experiment	
If not run for the full 64 hours, flow cells can be stored for subsequent flow cell checks and sequencing experiments. If the flow cell is ready for returning to Oxford Nanopore, please follow the returns and packing instructions. Flongle flow cells are single-use consumables. Once sequencing with a biological sample has occurred and the experiment has finished, dispose of the flow cell via a biological waste procedure.	