

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revU_14Aug2019
 Last update: 15/11/2023



Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 50 ng of poly(A)-tailed RNA or 500 ng of total RNA in 9 µl	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
<input type="checkbox"/> Direct RNA Sequencing Kit (SQK-RNA002)	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Magnetic rack, suitable for 1.5 ml Eppendorf tubes
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, AM9937)	<input type="checkbox"/> Microfuge
	<input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> SuperScript™ III Reverse Transcriptase (Thermo Fisher Scientific, cat # 18080044)	<input type="checkbox"/> Ice bucket with ice
	<input type="checkbox"/> 10 mM dNTP solution (e.g. NEB, cat # N0447)	<input type="checkbox"/> Timer
	<input type="checkbox"/> NEBNext® Quick Ligation Reaction Buffer (NEB, B6058)	<input type="checkbox"/> Thermal cycler
	<input type="checkbox"/> T4 DNA Ligase 2M U/ml (NEB, cat # M0202T/M)	<input type="checkbox"/> Qubit fluorometer (or equivalent for QC check)
	<input type="checkbox"/> Agencourt RNAClean XP beads (Beckman Coulter™, cat # A63987)	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> Qubit RNA HS Assay Kit (ThermoFisher, cat # Q32852)	
	<input type="checkbox"/> Qubit dsDNA HS Assay Kit (ThermoFisher, cat # Q32851)	
INSTRUCTIONS		NOTES/OBSERVATIONS
Library preparation		
Prepare the RNA in Nuclease-free water. <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 50 ng of poly(A)-tailed RNA or 500 ng of total RNA into a 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Adjust the volume to 9 µl with Nuclease-free water <input type="checkbox"/> Mix thoroughly by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge In a 0.2 ml thin-walled PCR tube, mix the reagents in the following order: <ul style="list-style-type: none"> <input type="checkbox"/> 3.0 µl NEBNext Quick Ligation Reaction Buffer <input type="checkbox"/> 9.0 µl RNA <input type="checkbox"/> 0.5 µl RNA CS (RCS), 110 nM <input type="checkbox"/> 1.0 µl RT Adapter (RTA) <input type="checkbox"/> 1.5 µl T4 DNA Ligase 		

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revU_14Aug2019
 Last update: 15/11/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Mix by pipetting and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</p> <p>Mix the following reagents together to make the reverse transcription master mix:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 9.0 µl Nuclease-free water <input type="checkbox"/> 2.0 µl 10 mM dNTPs <input type="checkbox"/> 8.0 µl 5x first-strand buffer <input type="checkbox"/> 4.0 µl 0.1 M DTT <p><input type="checkbox"/> Add the master mix to the 0.2 ml PCR tube containing the RT adapter-ligated RNA from the "RT Adapter ligation" step. Mix by pipetting.</p> <p><input type="checkbox"/> Add 2 µl of SuperScript III Reverse Transcriptase to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Place the tube in a thermal cycler and incubate at 50°C for 50 minutes, then 70°C for 10 minutes, and bring the sample to 4°C before proceeding to the next step.</p> <p><input type="checkbox"/> Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Resuspend the stock of Agencourt RNAClean XP beads by vortexing.</p> <p><input type="checkbox"/> Add 72 µl of resuspended Agencourt RNAClean XP beads to the reverse transcription reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 200 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant when clear and colourless.</p> <p>Keep the tube on magnet, and wash the beads with 150 µl of freshly prepared 70% ethanol without disturbing the pellet as described below.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Keeping the magnetic rack on the benchtop, rotate the bead-containing tube by 180°. Wait for the beads to migrate towards the magnet and form a pellet. <input type="checkbox"/> Rotate the tube 180° again (back to the starting position), and wait for the beads to pellet. <p><input type="checkbox"/> Remove the 70% ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet until the eluate is clear and colourless. Keep the tubes on the magnet and pipette off any residual ethanol.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 20 µl Nuclease-free water. Incubate for 5 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 20 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p>	

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revU_14Aug2019
 Last update: 15/11/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>In the same 1.5 ml Eppendorf DNA LoBind tube, mix the reagents in the following order:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 8.0 µl NEBNext Quick Ligation Reaction Buffer <input type="checkbox"/> 6.0 µl RNA Adapter (RMX) <input type="checkbox"/> 3.0 µl Nuclease-free water <input type="checkbox"/> 3.0 µl T4 DNA Ligase <p><input type="checkbox"/> Mix by pipetting.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</p> <p><input type="checkbox"/> Resuspend the stock of Agencourt RNAClean XP beads by vortexing.</p> <p><input type="checkbox"/> Add 16 µl of resuspended Agencourt RNAClean XP beads to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant when clear and colourless.</p> <p><input type="checkbox"/> Add 150 µl of the Wash Buffer (WSB) to the beads. Close the tube lid and resuspend the beads by flicking the tube. Return the tube to the magnetic rack, allow the beads to pellet for 5 minutes and pipette off the supernatant when clear and colourless.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 21 µl Elution Buffer by the gently flicking the tube. Incubate for 10 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 21 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Quantify 1 µl of reverse-transcribed and adapted RNA using the Qubit fluorometer DNA HS assay - recovery aim ~20 ng.</p>	
<p>The reverse-transcribed and adapted RNA is now ready for loading into the flow cell.</p>	
<p>Priming and loading the SpotON flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the RNA Running Buffer (RRB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT. <input type="checkbox"/> Mix the RNA Running Buffer (RRB), Flush Buffer (FB) and Flush Tether (FLT) tubes thoroughly by vortexing and spin down at RT. <input type="checkbox"/> To prepare the flow cell priming mix, add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing at RT. <input type="checkbox"/> Open the MinION device lid and slide the flow cell under the clip. <input type="checkbox"/> Slide the priming port cover clockwise to open the priming port. 	

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revU_14Aug2019
 Last update: 15/11/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>IMPORTANT</p> <p><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</p>	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <p><input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes. During this time, prepare the library for loading by following the steps below.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> Thoroughly mix the contents of the RRB tube by vortexing or pipetting, and spin down briefly.</p>	
<p><input type="checkbox"/> Take 20 µl of the prepared RNA library and mix it with 17.5 µl of Nuclease-free water.</p> <p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 37.5 µl RNA Running Buffer (RRB) <input type="checkbox"/> 37.5 µl RNA library in Nuclease-free water <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <p><input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading.</p> <p><input type="checkbox"/> Add 75 µl of sample to the Flow Cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</p> <p><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port and close the priming port.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output.</p>	
<p>Place the light shield onto the flow cell, as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Carefully place the leading edge of the light shield against the clip. Note: Do not force the light shield underneath the clip. <input type="checkbox"/> Gently lower the light shield onto the flow cell. The light shield should sit around the SpotON cover, covering the entire top section of the flow cell. 	
<p>Close the device lid and set up a sequencing run on MinKNOW.</p>	

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revU_14Aug2019
Last update: 15/11/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
Flow cell reuse and returns	
<input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C. <input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore.	
IMPORTANT <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.	