

Rapid sequencing gDNA - Field Sequencing Kit (SQK-LRK001)



Version: FSK\_9049\_v1\_revR\_14Aug2019  
 Last update: 07/09/2023

Flow Cell Number: .....

DNA Samples: .....

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> ~400 ng high molecular weight genomic DNA per sample		<input type="checkbox"/> Timer
<input type="checkbox"/> Field Sequencing Kit (SQK-LRK001)		<input type="checkbox"/> Method of achieving 80° C
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP001)		<input type="checkbox"/> Pipettes and pipette tips P10, P20, P100, P1000
INSTRUCTIONS		NOTES/OBSERVATIONS
<b>Library preparation</b>		
<p>DNA tagmentation</p> <p>Prepare the DNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Transfer ~400 ng genomic DNA into a DNA LoBind tube</li> <li><input type="checkbox"/> Adjust the volume to 10 µl with Nuclease-free water</li> <li><input type="checkbox"/> Mix by flicking the tube to avoid unwanted shearing</li> <li><input type="checkbox"/> Spin down briefly in a microfuge</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Cut off Tube 1 (left hand tube, containing FRL)</li> <li><input type="checkbox"/> Using a clean, empty pipette tip, pierce the foil of Tube 1. Take care to not disturb the pellet.</li> <li><input type="checkbox"/> Add 10 µl of the input DNA to Tube 1.</li> <li><input type="checkbox"/> Mix gently by pipetting up and down. Make sure all liquid is collected at the bottom of the tube.</li> <li><input type="checkbox"/> Incubate the tube at RT for 1 minute and then at 80° C for 1 minute.</li> </ul> <p>Adapter attachment</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Using a clean empty pipette tip, pierce the foil of Tube 2. Take care to not disturb the pellet.</li> <li><input type="checkbox"/> Transfer 10 µl of the tagmented DNA from Tube 1 to Tube 2.</li> <li><input type="checkbox"/> Mix gently by pipetting up and down. Make sure all liquid is collected at the bottom of the tube.</li> <li><input type="checkbox"/> Incubate the reaction for 5 minutes at RT.</li> <li><input type="checkbox"/> Using a clean empty pipette tip, pierce the foil of Tube 3 (right-hand tube). Take care to not disturb the pellet.</li> <li><input type="checkbox"/> Add 65 µl Resuspension Buffer (RTB) into Tube 3. Mix by pipetting up and down, and make sure all liquid is collected at the bottom of the tube.</li> </ul>		
The prepared DNA library is used for loading into the flow cell.		

Rapid sequencing gDNA - Field Sequencing Kit (SQK-LRK001)

Version: FSK\_9049\_v1\_revR\_14Aug2019  
 Last update: 07/09/2023



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Priming and loading the SpotON flow cell</b></p>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Mix the Flush Tether (FLT) and Flush Buffer (FLB) tubes by pipetting and spin down at RT.</li>   <li><input type="checkbox"/> Open the MinION Mk1B lid and slide the flow cell under the clip.</li>   <li><input type="checkbox"/> Slide the priming port cover clockwise to open the priming port.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 <math>\mu</math>l, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</li> </ul>	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 <math>\mu</math>l</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 <math>\mu</math>l, to draw back 20-30 <math>\mu</math>l, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Prepare the priming mix: add 30 <math>\mu</math>l of Flush Tether to the tube of Flush Buffer, and mix well by pipetting up and down.</li>   <li><input type="checkbox"/> Load 800 <math>\mu</math>l of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes. During this time, prepare the library for loading by following the steps below.</li>   <li><input type="checkbox"/> Prepare the library for loading: add 65 <math>\mu</math>l of the resuspended Tube 3 material into Tube 2. Mix by pipetting up and down.</li> </ul> <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li> <li><input type="checkbox"/> Load 200 <math>\mu</math>l of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li>   <li><input type="checkbox"/> Add 75 <math>\mu</math>l of the sample (mix from Tube 2 and Tube 3) to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</li>   <li><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION device lid.</li> </ul>	
<p><b>Flow cell reuse and returns</b></p>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C.</li> </ul>	

Rapid sequencing gDNA - Field Sequencing Kit (SQK-LRK001)



Version: FSK\_9049\_v1\_revR\_14Aug2019  
Last update: 07/09/2023

Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore.	
<b>IMPORTANT</b> <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.	