

Flow Cell Wash Kit (EXP-WSH003)

Version: WFC_9088_v1_revG_18Sep2019
 Last update: 04/07/2023



Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> Flow Cell Wash Kit (EXP-WSH003)		<input type="checkbox"/> Ice bucket with ice
		<input type="checkbox"/> Pipettes and pipette tips P20, P1000
INSTRUCTIONS		NOTES/OBSERVATIONS
Flushing a MinION/GridION Flow Cell		
<p>Preparation to run the washing procedure</p> <ul style="list-style-type: none"> <input type="checkbox"/> Place the tube of Wash Solution A on ice. Do not vortex the tube. <input type="checkbox"/> Thaw one tube of Wash Solution B at RT. <input type="checkbox"/> Mix the contents of Wash Solution B thoroughly by vortexing, spin down briefly and place on ice. <p>In a clean 1.5 ml Eppendorf DNA LoBind tube, prepare the following Wash Mix:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 20 µl Wash Solution A (A) <input type="checkbox"/> 380 µl Wash Solution B (B) <ul style="list-style-type: none"> <input type="checkbox"/> Mix well by pipetting, and place on ice. Do not vortex the tube. <input type="checkbox"/> Stop or pause the sequencing experiment in MinKNOW, and leave the flow cell in the device. <input type="checkbox"/> Before removing the waste fluid, ensure that the flow cell priming port cover and SpotON sample port cover are closed, as indicated in the figure below. <input type="checkbox"/> Remove all fluid from the waste channel through waste port 1 using a P1000 pipette. 		
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> It is vital that the flow cell priming port and SpotON sample port are closed before removing the waste buffer to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels. 		
<ul style="list-style-type: none"> <input type="checkbox"/> Slide the flow cell priming port cover clockwise to open. <p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl. <input type="checkbox"/> Insert the tip into the flow cell priming port. <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer/liquid entering the pipette tip. <input type="checkbox"/> Visually check that there is continuous buffer from the flow cell priming port across the sensor array. 		
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores. 		

Flow Cell Wash Kit (EXP-WSH003)

Version: WFC_9088_v1_revG_18Sep2019
 Last update: 04/07/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Load 400 µl of the prepared Flow Cell Wash Mix into the flow cell priming port, avoiding the introduction of air. <input type="checkbox"/> Close the priming port and wait for 30 minutes. <input type="checkbox"/> Before removing the waste fluid, ensure that the flow cell priming port cover and SpotON sample port cover are closed, as indicated in the figure below. <input type="checkbox"/> Remove all fluid from the waste channel through waste port 1 using a P1000 pipette. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> It is vital that the flow cell priming port and SpotON sample port are closed before removing the waste buffer to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels. 	
<p>Follow one of the two options described in the next steps of the protocol.</p>	
<p>To run a second library on a MinION/GridION flow cell straight away</p>	
<p>To run a second library straight away, follow the instructions in the "Priming and loading the flow cell" section of your library preparation protocol with the recommendations below.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Pipette very slowly when loading priming mix into the flow cell. <input type="checkbox"/> Wait five minutes between priming mix flushes. <input type="checkbox"/> After the five minute pause, close the priming port, ensure the SpotON port is closed and remove the waste from waste port 1. This prevents the nuclease from diffusing through the flow cell. Repeat this step after the second priming mix flush. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> When priming a flow cell after a nuclease wash with the Flow Cell Wash Kit, it is vital to wait five minutes between the priming mix flushes and to remove the waste for effective removal of the nuclease. 	
<p>To store the MinION/GridION flow cell for later use</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw one tube of Storage Buffer (S) at RT. <input type="checkbox"/> Mix contents thoroughly by pipetting and spin down briefly. <input type="checkbox"/> Slide the flow cell priming port cover clockwise to open. <p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl. <input type="checkbox"/> Insert the tip into the flow cell priming port. <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer/liquid entering the pipette tip. <input type="checkbox"/> Visually check that there is continuous buffer from the flow cell priming port across the sensor array. <ul style="list-style-type: none"> <input type="checkbox"/> Slowly add 500 µl of Storage Buffer (S) through the flow cell priming port. <input type="checkbox"/> Close the priming port. 	

Flow Cell Wash Kit (EXP-WSH003)

Version: WFC_9088_v1_revG_18Sep2019
Last update: 04/07/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Remove all fluid from the waste channel through waste port 1 using a P1000 pipette.	
IMPORTANT <input type="checkbox"/> It is vital that the flow cell priming port and SpotON sample port are closed before removing the waste buffer to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.	
<input type="checkbox"/> The flow cell can now be stored at 4-8°C.	
When you wish to reuse the flow cell, remove the flow cell from storage, and allow it to warm to RT for ~5 minutes.	