

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
 Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

### Before start checklist

Materials	Consumables	Equipment
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP004)	<input type="checkbox"/> MinION Flow Cell (FLO-MIN106 or FLO-MIN114)	<input type="checkbox"/> MinION or GridION device
<input type="checkbox"/> Sequencing Auxiliary Vials V14 (EXP-AUX003)	<input type="checkbox"/> Bovine Serum Albumin (BSA) (50 mg/ml) (e.g. Invitrogen™ UltraPure™ BSA 50 mg/ml, AM2616)	<input type="checkbox"/> P200 pipette and tips
<input type="checkbox"/> Flow Cell Wash Kit (EXP-WSH004) or Flow Cell Wash Kit XL (EXP-WSH004-XL)	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Ice bucket with ice
<input type="checkbox"/> SFB Expansion (EXP-SFB001)	<input type="checkbox"/> Agencourt AMPure XP beads (Beckman Coulter™ cat # A63881)	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> Qubit™ Assay Tubes (Invitrogen, Q32856)	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
	<input type="checkbox"/> Qubit dsDNA HS Assay Kit (Invitrogen, Q32851)	<input type="checkbox"/> Microfuge
		<input type="checkbox"/> Magnetic rack
		<input type="checkbox"/> Heating block
		<input type="checkbox"/> Qubit fluorometer (or equivalent)
		<input type="checkbox"/> Pipettes and pipette tips P20, P1000

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Transfer a library between MinION flow cells</b></p> <p>Preparation to transfer a library to a second flow cell</p> <p><b>IMPORTANT</b></p> <p><input type="checkbox"/> We recommend keeping the light shield on the original flow cell when recovering the library.</p> <p><input type="checkbox"/> Stop the sequencing run for the original flow cell on MinKNOW by clicking 'Stop'.</p> <p><input type="checkbox"/> Thaw and prepare the flow cell priming mix according to the "Priming and loading the SpotON flow cell" section of the suitable protocol.</p> <p><input type="checkbox"/> Open the MinION or GridION device lid and slide the second flow cell under the clip. Press down firmly on the flow cell to ensure correct thermal and electrical contact.</p> <p><input type="checkbox"/> To prime the second flow cell, slide the priming port cover clockwise to open the priming port.</p> <p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</p>	

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
 Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 µl</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Load 800 µl of the priming mix into the second flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes.</li> </ul> <p>Complete the flow cell priming for the second flow cell:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li> <li><input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li> <li><input type="checkbox"/> To prepare the original flow cell for library recovery, slide open the priming port cover and lift open the SpotON sample port cover.</li> <li><input type="checkbox"/> Set a pipette to 75 µl and fully depress the pipette before inserting the tip into the SpotON port of the original flow cell. Slowly aspirate to recover the DNA library from the flow cell.</li> <li><input type="checkbox"/> Add the recovered DNA library to the second flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</li> <li><input type="checkbox"/> Gently replace the SpotON sample port cover of the second flow cell, making sure the bung enters the SpotON port and close the priming port.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output.</li> </ul>	
<p>Place the light shield onto the flow cell, as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Carefully place the leading edge of the light shield against the clip. Note: Do not force the light shield underneath the clip.</li> <li><input type="checkbox"/> Gently lower the light shield onto the flow cell. The light shield should sit around the SpotON cover, covering the entire top section of the flow cell.</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> The original flow cell can be flushed with deionised water and returned to Oxford Nanopore.</li> <li><input type="checkbox"/> Start a new sequencing run on MinKNOW for the second flow cell.</li> </ul>	
<p>Using the suitable protocol for your DNA library, continue with the "Sequencing and data analysis" section to complete the experiment.</p>	
<p><b>Clean up and transfer a library between MinION Flow Cells</b></p>	
<p>Preparation to clean up a library before transfer to a second flow cell</p>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> We recommend keeping the light shield on the original flow cell when recovering the library.</li> </ul>	

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
 Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Thaw the kit components at RT and prepare as indicated by the table below:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> AMPure XP Beads (AXP)                             <ul style="list-style-type: none"> <li>- 1. Thaw at RT: ☑</li> <li>- 2. Mix well by vortexing: ☑</li> <li>- 3. Briefly spin down: X</li> <li>- 4. Keep on ice: X Keep at RT</li> </ul> </li> <li><input type="checkbox"/> Short Fragment Buffer (SFB)                             <ul style="list-style-type: none"> <li>- 1. Thaw at RT: ☑</li> <li>- 2. Mix well by vortexing: ☑</li> <li>- 3. Briefly spin down: ☑</li> <li>- 4. Keep on ice: ☑</li> </ul> </li> <li><input type="checkbox"/> Elution Buffer (EB)                             <ul style="list-style-type: none"> <li>- 1. Thaw at RT: ☑</li> <li>- 2. Mix well by vortexing: ☑</li> <li>- 3. Briefly spin down: ☑</li> <li>- 4. Keep on ice: ☑</li> </ul> </li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Stop the sequencing run for the original flow cell on MinKNOW by clicking 'Stop'.</li> <li><input type="checkbox"/> To prepare the original flow cell for library recovery, slide open the priming port cover and lift open the SpotON sample port cover.</li> <li><input type="checkbox"/> Set a pipette to 75 µl and fully depress the pipette before inserting the tip into the SpotON port of the original flow cell. Slowly aspirate to recover the DNA library from the flow cell.</li> <li><input type="checkbox"/> Transfer the recovered library to a fresh 1.5 ml Eppendorf DNA LoBind tube and store on ice.</li> <li><input type="checkbox"/> The original flow cell can be removed from the MinION or GridION device by sliding the flow cell from under the clip.</li> <li><input type="checkbox"/> The original flow cell can be flushed with deionised water and returned to Oxford Nanopore.</li> <li><input type="checkbox"/> Resuspend the AMPure XP Beads (AXP) by vortexing.</li> <li><input type="checkbox"/> Add 300 µl of resuspended AMPure XP Beads (AXP) to the recovered library and mix by flicking.</li> <li><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</li> <li><input type="checkbox"/> Spin down the sample and pellet on a magnet until supernatant is clear and colourless. Keep the tube on the magnet, and pipette off the supernatant.</li> <li><input type="checkbox"/> Wash the beads by adding 150 µl of Short Fragment Buffer (SFB). Flick the beads to resuspend, spin down, then return the tube to the magnetic rack and allow the beads to pellet. Remove the supernatant using a pipette and discard.</li> <li><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant.</li> <li><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 13 µl of Elution Buffer (EB).</li> <li><input type="checkbox"/> Spin down and incubate for 10 minutes at RT. For high molecular weight DNA, incubating at 37°C can improve recovery of long fragments.</li> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute.</li> <li><input type="checkbox"/> Remove and retain 13 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube.</li> </ul>	

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
 Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer. If the recovered library is below the detection level of the Qubit dsDNA HS Assay, we do not recommend continuing to load the flow cell.</p>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> The library can be stored at 4°C.</li> <li><input type="checkbox"/> Thaw and prepare the flow cell priming mix according to the "Priming and loading the SpotON flow cell" section of the suitable protocol.</li> <li><input type="checkbox"/> Open the MinION or GridION device lid and slide the second flow cell under the clip. Press down firmly on the flow cell to ensure correct thermal and electrical contact.</li> <li><input type="checkbox"/> To prime the second flow cell, slide the priming port cover clockwise to open the priming port.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</li> </ul>	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 µl</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Load 800 µl of the priming mix into the second flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes.</li> <li><input type="checkbox"/> Thoroughly mix the contents of the Library Beads/Loading Beads by pipetting.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> The Library Beads/Loading Beads tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use.</li> </ul>	
<p>In a new tube, prepare the recovered library for loading according to the "Priming and loading the SpotON flow cell" section of the suitable protocol to ensure you are using the correct reagents and volumes.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 37.5 µl Sequencing Buffer (SB)</li> <li><input type="checkbox"/> 25.5 µl Library Beads (LIB) mixed immediately before use, or Library Solution (LIS), if using</li> <li><input type="checkbox"/> 12 µl Recovered DNA library</li> </ul> <p>Complete the flow cell priming for the second flow cell:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li> <li><input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li> <li><input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading.</li> <li><input type="checkbox"/> Add the recovered DNA library to the second flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</li> </ul>	

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
 Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Gently replace the SpotON sample port cover of the second flow cell, making sure the bung enters the SpotON port and close the priming port.	
<p><b>IMPORTANT</b></p> <input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output.	
<p>Place the light shield onto the flow cell, as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Carefully place the leading edge of the light shield against the clip. Note: Do not force the light shield underneath the clip.</li> <li><input type="checkbox"/> Gently lower the light shield onto the flow cell. The light shield should sit around the SpotON cover, covering the entire top section of the flow cell.</li> </ul> <p><input type="checkbox"/> Start a new sequencing run on MinKNOW for the second flow cell.</p>	
<p>Using the suitable protocol for your DNA library, continue with the "Sequencing and data analysis" section to complete the experiment.</p>	
<p><b>Recover a library to replace on a washed MinION flow cell</b></p>	
<p>Preparation to recover and wash a library to replace on the same flow cell</p>	
<p><b>IMPORTANT</b></p> <input type="checkbox"/> We recommend keeping the light shield on the flow cell during library recovery, washing and reloading for optimal sequencing output.	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Pause the sequencing run for the original flow cell on MinKNOW by clicking 'Pause'.</li> <li><input type="checkbox"/> Place the tube of Wash Mix (WMX) on ice. Do not vortex the tube.</li> <li><input type="checkbox"/> Thaw one tube of Wash Diluent (DIL) at RT and mix the contents of Wash Diluent (DIL) thoroughly by vortexing. Then spin down briefly and place on ice.</li> <li><input type="checkbox"/> To prepare the original flow cell for library recovery, slide open the priming port cover.</li> </ul>	
<p><b>IMPORTANT</b></p> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 µl.</li> <li><input type="checkbox"/> Insert the tip into the flow cell priming port.</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer/liquid entering the pipette tip.</li> <li><input type="checkbox"/> Visually check that there is continuous buffer from the flow cell priming port across the sensor array.</li> </ul>	
<p><b>IMPORTANT</b></p> <input type="checkbox"/> Be aware that the library is removed from the priming port as a larger volume and expect to see the fluid in the waste channel to move back.	

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
 Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Set a pipette to 150 µl and fully depress the pipette before inserting the tip into the priming port of the original flow cell. Slowly aspirate to recover the DNA library from the flow cell.</p> <p><input type="checkbox"/> Transfer the recovered library to a fresh 1.5 ml Eppendorf DNA LoBind tube and store on ice.</p> <p>In a clean 1.5 ml Eppendorf DNA LoBind tube, prepare the following Flow Cell Wash Mix:</p> <p><input type="checkbox"/> 2 µl Wash Mix (WMX)</p> <p><input type="checkbox"/> 398 µl Wash Diluent (DIL)</p> <p><input type="checkbox"/> Mix well by pipetting, and place on ice. Do not vortex the tube.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> It is vital that the flow cell priming port and SpotON sample port are closed before removing the waste buffer to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.</p>	
<p>Remove the waste buffer, as follows:</p> <p><input type="checkbox"/> Close the priming port and SpotON sample port cover, as indicated in the figure below.</p> <p><input type="checkbox"/> Insert a P1000 pipette into waste port 1 and remove the waste buffer.</p> <p><input type="checkbox"/> Rotate the flow cell priming port cover clockwise so that the priming port is visible.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</p>	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <p><input type="checkbox"/> Set a P1000 pipette to 200 µl.</p> <p><input type="checkbox"/> Insert the tip into the flow cell priming port.</p> <p><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer/liquid entering the pipette tip.</p> <p><input type="checkbox"/> Visually check that there is continuous buffer from the flow cell priming port across the sensor array.</p> <p>Slowly load 200 µl of the prepared flow cell wash mix into the priming port, as follows:</p> <p><input type="checkbox"/> Using a P1000 pipette, take 200 µl of the flow cell wash mix</p> <p><input type="checkbox"/> Insert the pipette tip into the priming port, ensuring there are no bubbles in the tip</p> <p><input type="checkbox"/> Slowly twist the pipette wheel down to load the flow cell (if possible with your pipette) or push down the plunger very slowly, leaving a small volume of buffer in the pipette tip.</p> <p><input type="checkbox"/> Set a timer for a 5 minute incubation.</p> <p>Once the 5 minute incubation is complete, carefully load the remaining 200 µl of the prepared flow cell wash mix into the priming port, as follows:</p> <p><input type="checkbox"/> Using a P1000 pipette, take the remaining 200 µl of the flow cell wash mix</p> <p><input type="checkbox"/> Insert the pipette tip into the priming port, ensuring there are no bubbles in the tip</p> <p><input type="checkbox"/> Slowly twist the pipette wheel down to load the flow cell (if possible with your pipette) or push down the plunger very slowly, leaving a small volume of buffer in the pipette tip.</p> <p><input type="checkbox"/> Close the priming port and wait for 1 hour.</p>	

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
 Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> It is vital that the flow cell priming port and SpotON sample port are closed before removing the waste buffer to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.</p>	
<p>Remove the waste buffer, as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Close the priming port and SpotON sample port cover, as indicated in the figure below.</li> <li><input type="checkbox"/> Insert a P1000 pipette into waste port 1 and remove the waste buffer.</li> </ul> <p><input type="checkbox"/> Thaw and prepare the flow cell priming mix according to the "Priming and loading the SpotON flow cell" section of the suitable protocol.</p> <p><input type="checkbox"/> Slide the flow cell priming port cover clockwise to open the priming port.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</p>	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 µl</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <p><input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes.</p> <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li> <li><input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li> </ul> <p><input type="checkbox"/> Bring the recovered library up to RT and load 150 µl to the flow cell via the SpotON port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</p> <p><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port and close the priming port.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> If the light shield was removed during the washing step, the light shield should be replaced on the flow cell as soon as library is loaded for optimal sequencing output.</p>	
<p>Place the light shield onto the flow cell, as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Carefully place the leading edge of the light shield against the clip. Note: Do not force the light shield underneath the clip.</li> <li><input type="checkbox"/> Gently lower the light shield onto the flow cell. The light shield should sit around the SpotON cover, covering the entire top section of the flow cell.</li> </ul>	

# Library recovery from flow cells

Version: LIR\_9178\_v1\_revJ\_11Jan2023  
Last update: 28/02/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Restart the sequencing run on MinkNOW.	
Using the suitable protocol for your DNA library, continue with the "Sequencing and data analysis" section to complete the experiment.	