

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 100 ng of sheared genomic DNA or amplicon DNA	<input type="checkbox"/> MinION and GridION Flow Cell	<input type="checkbox"/> MinION or GridION device
<input type="checkbox"/> Ligation Sequencing Kit V14 (SQK-LSK114)	<input type="checkbox"/> LongAmp Hot Start Taq 2X Master Mix (NEB, M0533)	<input type="checkbox"/> MinION and GridION Flow Cell Light Shield
<input type="checkbox"/> PCR Expansion (EXP-PCA001)	<input type="checkbox"/> NEBNext Ultra II End repair/dA-tailing Module (NEB, E7546)	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
	<input type="checkbox"/> Salt-T4® DNA Ligase (NEB, M0467)	<input type="checkbox"/> Magnetic rack, suitable for 1.5 ml Eppendorf tubes
	<input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (NEB, M0367)	<input type="checkbox"/> Microfuge
	<input type="checkbox"/> Agencourt AMPure XP beads (Beckman Coulter, A63881)	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Thermal cycler
	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Ice bucket with ice
	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, AM9937)	<input type="checkbox"/> Timer
	<input type="checkbox"/> Freshly prepared 80% ethanol in nuclease-free water	<input type="checkbox"/> Qubit fluorometer (or equivalent for QC check)
	<input type="checkbox"/> Bovine Serum Albumin (BSA) (50 mg/ml) (e.g. Invitrogen™ UltraPure™ BSA 50 mg/ml, AM2616)	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
INSTRUCTIONS		NOTES/OBSERVATIONS
End-prep		
Check your flow cell.		
Prepare the NEBNext Ultra II End Repair / dA-tailing Module reagents in accordance with manufacturer's instructions, and place on ice: <ul style="list-style-type: none"> <input type="checkbox"/> Thaw all reagents on ice. <input type="checkbox"/> Ensure the reagents are well mixed. Note: Do not vortex the Ultra II End Prep Enzyme Mix. <input type="checkbox"/> Always spin down tubes before opening for the first time each day. <input type="checkbox"/> The NEBNext Ultra II End Prep Reaction Buffer may contain a white precipitate. If this occurs, allow the mixture(s) to come to RT and pipette the buffer several times to break up the precipitate, followed by a quick vortex to mix. 		

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Prepare the DNA in Nuclease-free water</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 100 ng genomic DNA into a 0.2 ml thin-walled PCR tube <input type="checkbox"/> Adjust the volume to 48 µl with Nuclease-free water <input type="checkbox"/> Mix thoroughly by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <p>In a 0.2 ml thin-walled PCR tube, mix the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 48 µl 100 ng DNA <input type="checkbox"/> 7 µl Ultra II End-prep Reaction Buffer <input type="checkbox"/> 3 µl Ultra II End-prep Enzyme Mix <input type="checkbox"/> 2 µl Nuclease-free water <ul style="list-style-type: none"> <input type="checkbox"/> Mix by pipetting and briefly spin down. <input type="checkbox"/> Using a thermal cycler, incubate at 20°C for 5 minutes and 65°C for 5 minutes. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Transfer the DNA sample to a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by flicking the tube. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 80% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet until supernatant is clear and colourless. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the sample from the magnet and resuspend the pellet in 16 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute. <input type="checkbox"/> Remove and retain 16 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Take forward the repaired and end-prepped DNA into the next step. However, at this point it is also possible to store the sample at 4°C overnight.</p>	

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>PCR adapters ligation and amplification</p> <p>Add the reagents as follows in the order below:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 15 µl End-prepped DNA <input type="checkbox"/> 10 µl PCR Adapter (PCA) <input type="checkbox"/> 25 µl Blunt/TA Ligase Master Mix <ul style="list-style-type: none"> <input type="checkbox"/> Mix by pipetting and briefly spin down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Add 20 µl of resuspended AMPure XP beads for a 0.4X clean and mix by pipetting up and down ten times. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 80% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant when clear and colourless. <input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 25 µl Nuclease-free water. Spin down and incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <p>Remove and retain the eluate once it is clear and colourless. Transfer the eluted sample to a fresh 0.2 ml PCR tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Dispose of the pelleted beads 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> If using amplicon input, start your second round of PCR from this point. Ensure the sample has undergone a round of PCR with tailed primers. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Quantify 1 µl of adapted DNA using a Qubit fluorometer. <p>Prepare the PCA adapted DNA in Nuclease-free water:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 20 ng genomic DNA into a 0.2 ml thin-walled PCR tube. <input type="checkbox"/> Adjust the volume to 48 µl with Nuclease-free water. <input type="checkbox"/> Mix thoroughly by flicking the tube. <input type="checkbox"/> Spin down briefly in a microfuge. 	

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Set up the adapted DNA PCR as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 46 µl Nuclease-free water <input type="checkbox"/> 2 µl Primer Mix (PRM) <input type="checkbox"/> 2 µl 10 ng/µl adapter ligated template <input type="checkbox"/> 50 µl HotStart LongAmp Taq 2x master mix <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Amplify using the following cycling conditions:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Initial denaturation 3 mins @ 94 °C (1 cycle) <input type="checkbox"/> Denaturation 15 secs @ 94 °C (15-18 (a) cycles) <input type="checkbox"/> Annealing 15 secs @ 56 °C (15-18 (a) cycles) <input type="checkbox"/> Extension 6 mins (b) @ 65 °C (15-18 (a) cycles) <input type="checkbox"/> Final extension 10 mins @ 65 °C (1 cycle) <input type="checkbox"/> Hold @ 4 °C <ul style="list-style-type: none"> <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 80% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant when clear and colourless. <input type="checkbox"/> Keep the samples on the magnet and wash the beads with 200 µl of freshly prepared 80% ethanol without disturbing the pellets. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 52 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 52 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Quantify 1 µl of adapted DNA using a Qubit fluorometer. 	
<p>Take forward the sample into the next step. However, at this point, it is also possible to store the sample at 4°C overnight.</p>	

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>End-prep</p> <p>Prepare the NEBNext Ultra II End Repair / dA-tailing Module reagents in accordance with manufacturer's instructions, and place on ice:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw all reagents on ice. <input type="checkbox"/> Ensure the reagents are well mixed. Note: Do not vortex the Ultra II End Prep Enzyme Mix. <input type="checkbox"/> Always spin down tubes before opening for the first time each day. <input type="checkbox"/> The NEBNext Ultra II End Prep Reaction Buffer may contain a white precipitate. If this occurs, allow the mixture(s) to come to RT and pipette the buffer several times to break up the precipitate, followed by a quick vortex to mix. <p>Prepare the PCR amplified DNA in Nuclease-free water</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 1 ug genomic DNA into a 0.2 ml thin-walled PCR tube. <input type="checkbox"/> Adjust the volume to 50 µl with Nuclease-free water. <input type="checkbox"/> Mix thoroughly by flicking the tube to avoid unwanted shearing. <input type="checkbox"/> Spin down briefly in a microfuge. <p>In a 0.2 ml thin-walled PCR tube, mix the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 50 µl DNA <input type="checkbox"/> 7 µl Ultra II End-prep Reaction Buffer <input type="checkbox"/> 3 µl Ultra II End-prep Enzyme Mix <ul style="list-style-type: none"> <input type="checkbox"/> Thoroughly mix the reaction by gently pipetting and briefly spinning down. <input type="checkbox"/> Using a thermal cycler, incubate at 20°C for 5 minutes and 65°C for 5 minutes. <input type="checkbox"/> Resuspend the AMPure XP Beads (AXP) by vortexing. <input type="checkbox"/> Transfer the DNA sample to a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Add 60 µl of resuspended the AMPure XP Beads (AXP) to the end-prep reaction and mix by flicking the tube. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 80% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet until supernatant is clear and colourless. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 61 µl Nuclease-free water. Incubate for 2 minutes at RT. 	

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute. <input type="checkbox"/> Remove and retain 61 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.	
Quantify 1 µl of eluted sample using a Qubit fluorometer.	
Take forward the repaired and end-prepped DNA into the adapter ligation step. However, at this point it is also possible to store the sample at 4°C overnight.	
Adapter ligation and clean-up	
<p>IMPORTANT</p> <input type="checkbox"/> Although third-party ligase products may be supplied with their own buffer, the ligation efficiency of the Ligation Adapter (LA) is higher when using the Ligation Buffer (LNB) supplied in the Ligation Sequencing Kit.	
<input type="checkbox"/> Spin down the Ligation Adapter (LA) and Salt-T4® DNA Ligase, and place on ice. <input type="checkbox"/> Thaw Ligation Buffer (LNB) at RT, spin down and mix by pipetting. Due to viscosity, vortexing this buffer is ineffective. Place on ice immediately after thawing and mixing. <input type="checkbox"/> Thaw the Elution Buffer (EB) at RT and mix by vortexing. Then spin down and place on ice.	
<p>IMPORTANT</p> Depending on the wash buffer (LFB or SFB) used, the clean-up step after adapter ligation is designed to either enrich for DNA fragments of >3 kb, or purify all fragments equally. <input type="checkbox"/> To enrich for DNA fragments of 3 kb or longer, use Long Fragment Buffer (LFB) <input type="checkbox"/> To retain DNA fragments of all sizes, use Short Fragment Buffer (SFB)	
<input type="checkbox"/> Thaw either Long Fragment Buffer (LFB) or Short Fragment Buffer (SFB) at RT and mix by vortexing. Then spin down and place on ice. In a 1.5 ml Eppendorf DNA LoBind tube, mix in the following order: <input type="checkbox"/> 60 µl DNA sample from the previous step <input type="checkbox"/> 5 µl Ligation Adapter (LA) <input type="checkbox"/> 25 µl Ligation Buffer (LNB) <input type="checkbox"/> 10 µl Salt-T4® DNA Ligase <input type="checkbox"/> Thoroughly mix the reaction by gently pipetting and briefly spinning down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. <input type="checkbox"/> Resuspend the AMPure XP Beads (AXP) by vortexing. <input type="checkbox"/> Add 40 µl of resuspended AMPure XP Beads (AXP) to the reaction and mix by flicking the tube. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant when clear and colourless.	

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Wash the beads by adding either 250 µl Long Fragment Buffer (LFB) or 250 µl Short Fragment Buffer (SFB). Flick the beads to resuspend, spin down, then return the tube to the magnetic rack and allow the beads to pellet. Remove the supernatant using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 15 µl Elution Buffer (EB). Spin down and incubate for 10 minutes at RT. For high molecular weight DNA, incubating at 37°C can improve the recovery of long fragments. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute. <input type="checkbox"/> Remove and retain 15 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	
<p>Depending on your DNA library fragment size, prepare your final library in 12 µl of Elution Buffer (EB).</p> <ul style="list-style-type: none"> <input type="checkbox"/> 100 fmol Very short (<1 kb) <input type="checkbox"/> 35–50 fmol Short (1–10 kb) <input type="checkbox"/> 300 ng Long (>10 kb) 	
<p>The prepared library is used for loading into the flow cell. Store the library on ice or at 4°C until ready to load.</p>	
<p>Priming and loading the MinION and GridION Flow Cell</p>	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Please note, this kit is only compatible with R10.4.1 flow cells (FLO-MIN114). 	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer (SB), Library Beads (LIB) or Library Solution (LIS, if using), Flow Cell Tether (FCT) and Flow Cell Flush (FCF) at RT before mixing by vortexing. Then spin down and store on ice. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> For optimal sequencing performance and improved output on MinION R10.4.1 flow cells (FLO-MIN114), we recommend adding Bovine Serum Albumin (BSA) to the flow cell priming mix at a final concentration of 0.2 mg/ml. 	
<p>To prepare the flow cell priming mix with BSA, combine Flow Cell Flush (FCF) and Flow Cell Tether (FCT), as directed below. Mix by pipetting at RT.</p> <ul style="list-style-type: none"> <input type="checkbox"/> 1,170 µl Flow Cell Flush (FCF) <input type="checkbox"/> 5 µl Bovine Serum Albumin (BSA) at 50 mg/ml <input type="checkbox"/> 30 µl Flow Cell Tether (FCT) <ul style="list-style-type: none"> <input type="checkbox"/> Open the MinION or GridION device lid and slide the flow cell under the clip. Press down firmly on the flow cell to ensure correct thermal and electrical contact. 	

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)

Version: LWP_9183_v114_revK_07Mar2023
Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Slide the flow cell priming port cover clockwise to open the priming port.	
<p>IMPORTANT</p> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes. During this time, prepare the library for loading by following the steps below. <input type="checkbox"/> Thoroughly mix the contents of the Library Beads (LIB) by pipetting. 	
<p>IMPORTANT</p> <input type="checkbox"/> The Library Beads (LIB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use.	
<p>In a new 1.5 ml Eppendorf DNA LoBind tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 37.5 µl Sequencing Buffer (SB) <input type="checkbox"/> 25.5 µl Library Beads (LIB) mixed immediately before use, or Library Solution (LIS), if using <input type="checkbox"/> 12 µl DNA library <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of the prepared library to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port and close the priming port. 	
<p>IMPORTANT</p> <input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output.	
<p>Place the light shield onto the flow cell, as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Carefully place the leading edge of the light shield against the clip. Note: Do not force the light shield underneath the clip. <input type="checkbox"/> Gently lower the light shield onto the flow cell. The light shield should sit around the SpotON cover, covering the entire top section of the flow cell. 	

Ligation sequencing V14 - low input by PCR (SQK-LSK114 with EXP-PCA001)



Version: LWP_9183_v114_revK_07Mar2023
 Last update: 21/08/2024

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
Close the device lid and set up a sequencing run on MinKNOW.	
Flow cell reuse and returns	
<input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C. <input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore.	
IMPORTANT <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.	