

Ligation sequencing gDNA - Multiplex Ligation Sequencing Kit  
V14 XL (SQK-MLK114.96-XL)



Version: MLK\_9193\_v114\_revA\_20Sep2023  
Last update: 04/10/2023

Flow Cell Number: .....

DNA Samples: .....

**Before start checklist**

**Materials**

- Multiplex Ligation Sequencing Kit V14 XL (SQK-MLK114.96-XL)
- 1000 ng gDNA per sample

**Consumables**

- NEB Blunt/TA Ligase Master Mix (NEB, M0367)
- NEBNext FFPE Repair Mix (NEB, M6630)
- NEBNext Ultra II End repair/dA-tailing Module (NEB, E7546)
- NEBNext Quick Ligation Module (NEB, E6056)
- Eppendorf twin.tec® PCR plate 96 LoBind, semi-skirted (Eppendorf™, cat # 0030129504) with heat seals
- 0.2 ml thin-walled PCR tubes
- 1.5 ml Eppendorf DNA LoBind tubes
- 2 ml Eppendorf DNA LoBind tubes
- Nuclease-free water (e.g. ThermoFisher, AM9937)
- Freshly prepared 80% ethanol in nuclease-free water
- Qubit™ Assay Tubes (Invitrogen, Q32856)
- Qubit dsDNA HS Assay Kit (Invitrogen, Q32851)

**Equipment**

- Hula mixer (gentle rotator mixer)
- Microfuge
- Magnetic rack
- Magnetic rack
- Vortex mixer
- Thermal cycler
- Ice bucket with ice
- Timer
- Qubit fluorometer (or equivalent)
- Pipettes and pipette tips P2, P10, P20, P100, P200, P1000, Multichannel

**INSTRUCTIONS**

**NOTES/OBSERVATIONS**

**DNA repair and end-prep**

- Thaw the AMPure XP Beads (AXP) at RT and mix by vortexing. Keep the beads at RT.
- Prepare the NEBNext FFPE DNA Repair Mix and NEBNext Ultra II End Repair / dA-tailing Module reagents in accordance with manufacturer's instructions, and place on ice.
- Thaw all reagents on ice.
  - Flick and/or invert the reagent tubes to ensure they are well mixed.  
Note: Do not vortex the FFPE DNA Repair Mix or Ultra II End Prep Enzyme Mix.
  - Always spin down tubes before opening for the first time each day.
  - The Ultra II End Prep Buffer and FFPE DNA Repair Buffer may have a little precipitate. Allow the mixture to come to RT and pipette the buffer up and down several times to break up the precipitate, followed by vortexing the tube for 30 seconds to solubilise any precipitate.  
Note: It is important the buffers are mixed well by vortexing.
  - The FFPE DNA Repair Buffer may have a yellow tinge and is fine to use if yellow.

Ligation sequencing gDNA - Multiplex Ligation Sequencing Kit  
V14 XL (SQK-MLK114.96-XL)



Version: MLK\_9193\_v114\_revA\_20Sep2023  
Last update: 04/10/2023

Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Do not vortex the NEBNext FFPE DNA Repair Mix or NEBNext Ultra II End Prep Enzyme Mix.</p> <p><input type="checkbox"/> In clean 0.2 ml thin-walled PCR tubes, aliquot 1000 ng per sample.</p> <p><input type="checkbox"/> Make up each sample to 12 µl using Nuclease-free water. Mix gently by pipetting and spin down.</p> <p>Combine the following components per sample:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 12 µl DNA sample</li> <li><input type="checkbox"/> 0.875 µl NEBNext FFPE DNA Repair Buffer</li> <li><input type="checkbox"/> 0.875 µl Ultra II End-prep reaction buffer</li> <li><input type="checkbox"/> 0.75 µl Ultra II End-prep enzyme mix</li> <li><input type="checkbox"/> 0.50 µl NEBNext FFPE DNA Repair Mix</li> </ul> <p><input type="checkbox"/> Ensure the components are thoroughly mixed by pipetting and spin down briefly.</p> <p><input type="checkbox"/> Using a thermal cycler, incubate at 20°C for 5 minutes and 65°C for 5 minutes.</p> <p><input type="checkbox"/> Transfer each sample to clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Resuspend the AMPure XP beads (AXP) by vortexing.</p> <p><input type="checkbox"/> Add 15 µl of resuspended AMPure XP Beads (AXP) to each end-prep reaction and mix by flicking the tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of 80% ethanol in Nuclease-free water per sample.</p> <p><input type="checkbox"/> Spin down the samples and pellet the beads on a magnet until the eluate is clear and colourless. Keep the tubes on the magnet and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tubes on the magnet and wash the beads with 200 µl of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Briefly spin down and place the tubes back on the magnet. Pipette off any residual ethanol. Allow to dry for 30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tubes from the magnetic rack and resuspend the pellet in 10 µl Nuclease-free water. Spin down and incubate for 2 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 10 µl of eluate for each sample into clean 1.5 ml Eppendorf DNA LoBind tubes, individually.</p>	
<p>Quantify 1 µl of each eluted sample using a Qubit fluorometer.</p>	
<p><input type="checkbox"/> Keeping your samples separate, standardise them to an equimolar mass.</p>	

Ligation sequencing gDNA - Multiplex Ligation Sequencing Kit  
V14 XL (SQK-MLK114.96-XL)

Version: MLK\_9193\_v114\_revA\_20Sep2023  
Last update: 04/10/2023



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Take forward the equimolar samples in 7.5 µl to be barcoded and pooled in the native barcode ligation step. However, at this point it is also possible to store the sample at 4°C overnight.</p>	
<p><b>Native barcode ligation</b></p>	
<p>Prepare the NEB Blunt/TA Ligase Master Mix according to the manufacturer's instructions, and place on ice:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Thaw the reagents at RT.</li> <li><input type="checkbox"/> Spin down the reagent tubes for 5 seconds.</li> <li><input type="checkbox"/> Ensure the reagents are fully mixed by performing 10 full volume pipette mixes.</li> </ul> <p><input type="checkbox"/> Thaw the EDTA at RT and mix by vortexing. Then spin down and place on ice.</p> <p><input type="checkbox"/> Thaw the native barcodes at RT. Use one barcode per sample. Individually mix the barcodes by pipetting, spin down, and place them on ice.</p> <p>Select a unique barcode for each sample to be run in a group:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> For two samples on one flow cell, select two unique barcodes, one for each sample.</li> <li><input type="checkbox"/> For three samples on two flow cells, select three unique barcodes, one for each sample.</li> </ul> <p>In clean 1.5 ml Eppendorf DNA LoBind tubes, add the reagents in the following order per sample:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 7.5 µl End-prepped DNA</li> <li><input type="checkbox"/> 2.5 µl Native barcode</li> <li><input type="checkbox"/> 10 µl Blunt/TA Ligase Master Mix</li> </ul> <p><input type="checkbox"/> Thoroughly mix the reaction by gently pipetting and briefly spinning down.</p> <p><input type="checkbox"/> Incubate for 20 minutes at RT.</p> <p><input type="checkbox"/> Add 4 µl of EDTA to each tube and mix thoroughly by pipetting and spin down briefly.</p> <p><input type="checkbox"/> Pool each group of two or three uniquely barcoded samples in a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Ensure the beads are at RT and resuspend the AMPure XP beads (AXP) by vortexing.</p> <p>Add AMPure XP Beads (AXP) to the pooled reaction in a 0.4X ratio and mix by pipetting. The volume for AMPure XP Beads (AXP) will vary depending on the number of barcoded samples in the pool:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> For two barcoded samples add 19 µl of AMPure XP Beads (AXP)</li> <li><input type="checkbox"/> For three barcoded samples add 29 µl of AMPure XP Beads (AXP)</li> </ul> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p> <p><input type="checkbox"/> Prepare 1 ml of fresh 80% ethanol in Nuclease-free water per barcoded sample pool.</p> <p><input type="checkbox"/> Spin down the sample for 5 seconds and pellet on a magnet for 5 minutes. Keep the tube on the magnetic rack until the eluate is clear and colourless, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tube on the magnetic rack and wash the beads with 200 µl of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p>	

Ligation sequencing gDNA - Multiplex Ligation Sequencing Kit  
V14 XL (SQK-MLK114.96-XL)



Version: MLK\_9193\_v114\_revA\_20Sep2023  
Last update: 04/10/2023

Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <li><input type="checkbox"/> Spin down and place the tube back on the magnetic rack. Pipette off any residual ethanol. Allow the pellet to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</li> <li><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 35 µl Nuclease-free water by gently flicking.</li> <li><input type="checkbox"/> Incubate for 10 minutes at 37°C. Every 2 minutes, agitate the sample by gently flicking for 10 seconds to encourage DNA elution.</li> <li><input type="checkbox"/> Spin down the tube for 5 seconds and pellet the beads on a magnetic rack until the eluate is clear and colourless.</li> <li><input type="checkbox"/> Remove and retain 35 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</li> </ul>	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	
<p>Take forward the barcoded DNA library to the adapter ligation and clean-up step. However, you may store the sample at 4°C overnight.</p>	
<p><b>Adapter ligation and clean-up</b></p>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> The Native Adapter (NA) used in this kit and protocol is not interchangeable with other sequencing adapters.</li> </ul>	
<p>Prepare the NEBNext Quick Ligation Reaction Module according to the manufacturer's instructions, and place on ice:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Thaw the reagents at RT.</li> <li><input type="checkbox"/> Spin down the reagent tubes for 5 seconds.</li> <li><input type="checkbox"/> Ensure the reagents are fully mixed by performing 10 full volume pipette mixes.</li> </ul> <p>Note: Do NOT vortex the Quick T4 DNA Ligase.</p>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Do not vortex the Quick T4 DNA Ligase.</li> </ul>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Spin down the Native Adapter (NA) and Quick T4 DNA Ligase, pipette mix and place on ice.</li> <li><input type="checkbox"/> Thaw the Elution Buffer (EB) at RT and mix by vortexing. Then spin down and place on ice.</li> <li><input type="checkbox"/> To enrich for DNA fragments of 3 kb or longer, thaw one tube of Long Fragment Buffer (LFB) at RT, mix by vortexing, spin down and place on ice.</li> </ul> <p>In a 1.5 ml Eppendorf LoBind tube, mix in the following order:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 30 µl Pooled barcoded sample</li> <li><input type="checkbox"/> 5 µl Native Adapter (NA)</li> <li><input type="checkbox"/> 10 µl NEBNext Quick Ligation Reaction Buffer (5X)</li> <li><input type="checkbox"/> 5 µl Quick T4 DNA Ligase</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Thoroughly mix the reaction by gently pipetting and briefly spinning down.</li> <li><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</li> </ul>	

Ligation sequencing gDNA - Multiplex Ligation Sequencing Kit  
V14 XL (SQK-MLK114.96-XL)

Version: MLK\_9193\_v114\_revA\_20Sep2023  
Last update: 04/10/2023



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> The next clean-up step uses Long Fragment Buffer (LFB) rather than 80% ethanol to wash the beads. The use of ethanol will be detrimental to the sequencing reaction.</p>	
<p><input type="checkbox"/> Resuspend the AMPure XP Beads (AXP) by vortexing.</p> <p><input type="checkbox"/> Add 20 µl of resuspended AMPure XP Beads (AXP) to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p> <p><input type="checkbox"/> Spin down the sample and pellet on the magnetic rack. Keep the tube on the magnet and pipette off the supernatant.</p> <p><input type="checkbox"/> Wash the beads by adding 125 µl Long Fragment Buffer (LFB). Flick the beads to resuspend, spin down, then return the tube to the magnetic rack and allow the beads to pellet. Remove the supernatant using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 35 µl Elution Buffer (EB).</p> <p><input type="checkbox"/> Spin down and incubate for 10 minutes at 37°C. Every 2 minutes, agitate the sample by gently flicking for 10 seconds to encourage DNA elution.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute.</p> <p><input type="checkbox"/> Remove and retain 35 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube.</p>	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	
<p>Make up each library to 32 µl at 50 fmol, using Elution Buffer (EB).</p> <p><input type="checkbox"/> For two samples on one flow cell, make one DNA library from your eluate.</p> <p><input type="checkbox"/> For three samples on two flow cells, make two DNA libraries from your eluate.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Where possible, we recommend loading ~50 fmol of this final prepared library onto the R10.4.1 flow cell for our Multiplex Ligation Sequencing V14 protocol.</p>	
<p>The prepared library is used for loading into the flow cell. Store the library on ice or at 4°C until ready to load.</p>	
<p><b>Priming and loading multiple flow cells on a PromethION</b></p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> This kit is only compatible with R10.4.1 flow cells (FLO-PRO114M).</p>	
<p><input type="checkbox"/> Thaw the Sequencing Buffer (SB), Library Beads (LIB) or Library Solution (LIS, if using), Flow Cell Tether (FCT) and Flow Cell Flush (FCF) at RT, before mixing by vortexing. Then spin down before storing on ice.</p>	

Ligation sequencing gDNA - Multiplex Ligation Sequencing Kit  
V14 XL (SQK-MLK114.96-XL)



Version: MLK\_9193\_v114\_revA\_20Sep2023  
Last update: 04/10/2023

Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Scale up reagent volumes as needed.</p>	
<p>Prepare the flow cell priming mix in a suitable tube for the number of flow cells to flush. Once combined, mix well by briefly vortexing.</p> <p><input type="checkbox"/> 1,170 µl Flow Cell Flush (FCF)</p> <p><input type="checkbox"/> 30 µl Flow Cell Tether (FCT)</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> After taking flow cells out of the fridge, wait 20 minutes before inserting the flow cell into the PromethION for the flow cell to come to RT. Condensation can form on the flow cell in humid environments. Inspect the gold connector pins on the top and underside of the flow cell for condensation and wipe off with a lint-free wipe if any is observed. Ensure the heat pad (black pad) is present on the underside of the flow cell.</p>	
<p>For PromethION 2 Solo, load the flow cell(s) as follows:</p> <p><input type="checkbox"/> Place the flow cell flat on the metal plate.</p> <p><input type="checkbox"/> Slide the flow cell into the docking port until the gold pins or green board cannot be seen.</p> <p>For the PromethION 24/48, load the flow cell(s) into the docking ports:</p> <p><input type="checkbox"/> Line up the flow cell with the connector horizontally and vertically before smoothly inserting into position.</p> <p><input type="checkbox"/> Press down firmly onto the flow cell and ensure the latch engages and clicks into place.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Insertion of the flow cells at the wrong angle can cause damage to the pins on the PromethION and affect your sequencing results. If you find the pins on a PromethION position are damaged, please contact support@nanoporetech.com for assistance.</p>	
<p><input type="checkbox"/> If not already completed, perform a flow cell check on all flow cells.</p> <p><input type="checkbox"/> Slide the inlet port cover clockwise to open.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</p>	
<p>After opening the inlet port, draw back a small volume to remove any air bubbles:</p> <p><input type="checkbox"/> Set a P1000 pipette tip to 200 µl.</p> <p><input type="checkbox"/> Insert the tip into the inlet port.</p> <p><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you see a small volume of buffer entering the pipette tip.</p> <p><input type="checkbox"/> Load 500 µl of the priming mix into the flow cell via the inlet port, avoiding the introduction of air bubbles. Wait five minutes. During this time, prepare the library for loading using the next steps in the protocol.</p> <p><input type="checkbox"/> Thoroughly mix the contents of the Library Beads (LIB) by pipetting.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> The Library Beads (LIB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use.</p>	

Ligation sequencing gDNA - Multiplex Ligation Sequencing Kit  
 V14 XL (SQK-MLK114.96-XL)



Version: MLK\_9193\_v114\_revA\_20Sep2023  
 Last update: 04/10/2023

Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>In a new 1.5 ml Eppendorf DNA LoBind tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 100 µl Sequencing Buffer (SB)</li> <li><input type="checkbox"/> 68 µl Library Beads (LIB) thoroughly mixed before use, or Library Solution (LIS)</li> <li><input type="checkbox"/> 32 µl DNA library</li> </ul> <p><input type="checkbox"/> Complete the flow cell priming by slowly loading 500 µl of the priming mix into the inlet port.</p> <p><input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading.</p> <p><input type="checkbox"/> Load 200 µl of library into the inlet port using a P1000 pipette.</p> <p><input type="checkbox"/> Close the valve to seal the inlet port and close the PromethION lid when ready.</p>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output.</li> </ul>	
<p>If the light shield has been removed from the flow cell, install the light shield as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Align the inlet port cut out of the light shield with the inlet port cover on the flow cell. The leading edge of the light shield should sit above the flow cell ID.</li> <li><input type="checkbox"/> Firmly press the light shield around the inlet port cover. The inlet port clip will click into place underneath the inlet port cover.</li> </ul> <p><input type="checkbox"/> For multiple flow cell washing, use the same experiment name and identifying sample IDs for all runs to enable all flow cells to be paused simultaneously.</p>	
<p><b>Flow cell reuse and returns</b></p>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C.</li> <li><input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.</li> </ul>	