

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> Input influenza RNA	<input type="checkbox"/> SuperScript III One-Step RT-PCR System with Platinum Taq DNA Polymerase (ThermoFisher, cat # 12574018 or 12574026)	<input type="checkbox"/> Magnetic rack, suitable for 1.5 ml Eppendorf tubes
<input type="checkbox"/> Influenza A primers	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, AM9937)	<input type="checkbox"/> Magnetic rack suitable for 96 well plates, e.g. DynaMag™-96 Side Skirted Magnet (Thermo Fisher CAT#12027)
<input type="checkbox"/> Influenza B primers	<input type="checkbox"/> Agencourt AMPure XP beads (Beckman Coulter™ cat # A63881)	<input type="checkbox"/> Microfuge
<input type="checkbox"/> Ligation Sequencing Kit (SQK-LSK109)	<input type="checkbox"/> Freshly prepared 80% ethanol in nuclease-free water	<input type="checkbox"/> Vortex mixer
<input type="checkbox"/> Native Barcoding Expansion 96 (EXP-NBD196)	<input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (NEB, M0367)	<input type="checkbox"/> Thermal cycler
<input type="checkbox"/> Or Native Barcoding Expansion 1-12 (EXP-NBD104) and 13-24 (EXP-NBD114) if multiplexing 1-24 samples.	<input type="checkbox"/> NEBNext Ultra II End repair/dA-tailing Module (NEB, E7546)	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> NEBNext Quick Ligation Module (NEB, E6056)	<input type="checkbox"/> Microplate centrifuge, e.g. Fisherbrand™ Mini Plate Spinner Centrifuge (Fisher Scientific, 11766427)
<input type="checkbox"/> Adapter Mix II Expansion (EXP-AMII001)	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Qubit fluorometer (or equivalent for QC check)
<input type="checkbox"/> SFB Expansion (EXP-SFB001)	<input type="checkbox"/> 5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Multichannel pipettes suitable for dispensing 0.5–10 µl, 2–20 µl and 20–200 µl, and tips
<input type="checkbox"/> Sequencing Auxiliary Vials (EXP-AUX001)	<input type="checkbox"/> 15 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Ice bucket with ice
	<input type="checkbox"/> Reagent reservoirs for multichannel pipetting	<input type="checkbox"/> Timer
	<input type="checkbox"/> Qubit dsDNA HS Assay Kit (Invitrogen, Q32851)	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> Qubit™ Assay Tubes (Invitrogen, Q32856)	
	<input type="checkbox"/> Eppendorf twin.tec® PCR plate 96 LoBind, semi-skirted (Cat # 0030129504) with heat seals	
INSTRUCTIONS		NOTES/OBSERVATIONS
Reverse transcription, PCR and clean-up		
IMPORTANT		
<input type="checkbox"/> Keep the RNA sample on ice as much as possible to prevent nucleolytic degradation, which may affect sensitivity.		

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>To reduce risk of contamination, we recommend the use of PCR hoods with a UV steriliser when setting up the PCR plates.</p> <p>In a clean template-free pre-PCR hood, prepare the primer mixes for influenza A and influenza B as follows in 1.5 ml Eppendorf DNA LoBind tubes:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Nuclease-free water <ul style="list-style-type: none"> - Concentration: - - Volume: 378 µl <input type="checkbox"/> Tuni 12 <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 16.8 µl <input type="checkbox"/> Tuni 12.4 <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 4.2 µl <input type="checkbox"/> Tuni 13 <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 21 µl <input type="checkbox"/> Nuclease-free water <ul style="list-style-type: none"> - Concentration: - - Volume: 378 µl <input type="checkbox"/> B-PBs-UniF <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 5 µl <input type="checkbox"/> B-PBs-UniR <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 5 µl <input type="checkbox"/> B-PA-UniF <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 2.5 µl <input type="checkbox"/> B-PA-UniR <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 2.5 µl <input type="checkbox"/> B-HANA-UniF <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 5 µl <input type="checkbox"/> B-HANA-UniR <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 5 µl <input type="checkbox"/> B-NP-UniF <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 3 µl <input type="checkbox"/> B-NP-UniR <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 3 µl <input type="checkbox"/> B-M-Uni3F <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 1.5 µl <input type="checkbox"/> B-Mg-Uni3F <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 1.5 µl <input type="checkbox"/> B-M-Uni3R <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 3 µl <input type="checkbox"/> B-NS-Uni3F <ul style="list-style-type: none"> - Concentration: 100 µM - Volume: 2.5 µl 	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>In the template-free pre-PCR hood, prepare the following master mixes in Eppendorf DNA LoBind tubes and mix thoroughly as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Nuclease free water <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 280 µl - Influenza B RT-PCR Master Mix: 280 µl <input type="checkbox"/> Influenza A primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 28 µl - Influenza B RT-PCR Master Mix: - <input type="checkbox"/> Influenza B primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: - - Influenza B RT-PCR Master Mix: 28 µl <input type="checkbox"/> 2X Reaction Mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 350 µl - Influenza B RT-PCR Master Mix: 350 µl <input type="checkbox"/> SuperScript™ III RT/Platinum™ Taq Mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 28 µl - Influenza B RT-PCR Master Mix: 28 µl <input type="checkbox"/> Nuclease free water <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 560 µl - Influenza B RT-PCR Master Mix: 560 µl <input type="checkbox"/> Influenza A primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 56 µl - Influenza B RT-PCR Master Mix: - <input type="checkbox"/> Influenza B primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: - - Influenza B RT-PCR Master Mix: 56 µl <input type="checkbox"/> 2X Reaction Mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 700 µl - Influenza B RT-PCR Master Mix: 700 µl <input type="checkbox"/> SuperScript™ III RT/Platinum™ Taq Mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 56 µl - Influenza B RT-PCR Master Mix: 56 µl <input type="checkbox"/> Nuclease free water <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 1120 µl - Influenza B RT-PCR Master Mix: 1120 µl <input type="checkbox"/> Influenza A primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 112 µl - Influenza B RT-PCR Master Mix: - <input type="checkbox"/> Influenza B primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: - - Influenza B RT-PCR Master Mix: 112 µl <input type="checkbox"/> 2X Reaction Mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 1400 µl - Influenza B RT-PCR Master Mix: 1400 µl <input type="checkbox"/> SuperScript™ III RT/Platinum™ Taq Mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 112 µl - Influenza B RT-PCR Master Mix: 112 µl <input type="checkbox"/> Nuclease free water <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 2240 µl - Influenza B RT-PCR Master Mix: 2240 µl <input type="checkbox"/> Influenza A primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: 224 µl - Influenza B RT-PCR Master Mix: - <input type="checkbox"/> Influenza B primer mix <ul style="list-style-type: none"> - Influenza A RT-PCR Master Mix: - - Influenza B RT-PCR Master Mix: 224 µl 	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> For each influenza type, place one clean 96-well RT-PCR plate into a PCR-cooler (if using). <input type="checkbox"/> Using a stepper pipette or a multichannel pipette, aliquot 49 µl of influenza A RT-PCR Master Mix into the influenza A RT-PCR plate. <input type="checkbox"/> Using a stepper pipette or a multichannel pipette, aliquot 49 µl of influenza B RT-PCR Master Mix into the influenza B RT-PCR plate. <input type="checkbox"/> Seal the RT-PCR plate(s) and transfer to a template-addition pre-PCR hood. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> We recommend having a negative control for every plate of samples. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Transfer 1 µl of influenza A samples to the wells containing influenza A RT-PCR Master Mix in the influenza A RT-PCR plate and mix thoroughly by pipetting the contents of each well up and down. <input type="checkbox"/> Transfer 1 µl of influenza B samples to the wells containing influenza B RT-PCR Master Mix in the influenza B RT-PCR plate and mix thoroughly by pipetting the contents of each well up and down. <input type="checkbox"/> Seal the RT-PCR plate(s) and spin down in a centrifuge. <p>Incubate the influenza A RT-PCR plate using the following program, with the heated lid set to 105°C:</p> <ul style="list-style-type: none"> <input type="checkbox"/> cDNA synthesis 60 min @ 42°C (1 cycle) <input type="checkbox"/> Initial denaturation 2 min @ 94°C (1 cycle) <input type="checkbox"/> Denaturation Annealing and extension 30 sec 30 sec 3 min @ 94°C 45°C 68°C (5 cycles) <input type="checkbox"/> Denaturation Annealing and extension 30 sec 30 sec 3 min @ 94°C 57°C 68°C (31 cycles) <input type="checkbox"/> Hold @ 4°C <p>Incubate the influenza B RT-PCR plate using the following program, with the heated lid set to 105°C:</p> <ul style="list-style-type: none"> <input type="checkbox"/> cDNA synthesis 60 min @ 45°C (1 cycle) <input type="checkbox"/> cDNA synthesis 30 min @ 55°C (1 cycle) <input type="checkbox"/> Initial denaturation 2 min @ 94°C (1 cycle) <input type="checkbox"/> Denaturation Annealing and extension 20 sec 30 sec 3 min 30 sec @ 94°C 40°C 68°C (5 cycles) <input type="checkbox"/> Denaturation Annealing and extension 20 sec 30 sec 3 min 30 sec @ 94°C 58°C 68°C (30 cycles) <input type="checkbox"/> Final extension 10 min @ 68°C (1 cycle) <input type="checkbox"/> Hold @ 4°C 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> If available, a clean post-PCR hood should be used for all steps that involve handling amplified material. Decontamination with UV and or DNAzap between sample batches is recommended. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Add 50 µl of resuspended AMPure XP beads to each well of the RT-PCR plate(s) and mix by gently pipetting. <input type="checkbox"/> Incubate the PT-PCR plate(s) at RT for 10 minutes. <input type="checkbox"/> Prepare at least 500 µl 80% ethanol in Nuclease-free water per sample. 	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Spin down the RT-PCR plate(s) and pellet the beads on a magnet for 5 minutes. Keep the plate on the magnet until the eluate is clear and colourless, and pipette off the supernatant. <input type="checkbox"/> Keep the plate on the magnet and wash the beads in each well with 200 µl of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the plate back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the plate from the magnetic rack and resuspend each pellet in 15 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 15 µl of eluate containing the DNA per well, into a clean 96-well plate(s). 	
<p>Quantify 1 µl of each eluted sample using a Qubit fluorometer.</p>	
<p>Take forward your quantified sample to the end-prep step.</p>	
<p>End-prep</p>	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> We recommended carrying the negative control through this step until sequencing. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Determine the volume of the cleaned-up PCR reaction that yields 200 fmol of DNA per sample and aliquot into a clean 96-well plate (End-prep plate). <p>Prepare the NEBNext Ultra II End Repair / dA-tailing Module reagents in accordance with manufacturer's instructions, and place on ice:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw all reagents on ice. <input type="checkbox"/> Flick and/or invert the reagent tubes to ensure they are well mixed. Note: Do not vortex the Ultra II End Prep Enzyme Mix. <input type="checkbox"/> Always spin down tubes before opening for the first time each day. <input type="checkbox"/> The Ultra II End Prep Buffer may have a little precipitate. Allow the mixture to come to RT and pipette the buffer up and down several times to break up the precipitate, followed by vortexing the tube for 30 seconds to solubilise any precipitate. <ul style="list-style-type: none"> <input type="checkbox"/> Make up each sample per well to 12.5 µl using Nuclease-free water. 	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Prepare the following end-prep master mix in 1.5 ml Eppendorf DNA LoBind tube and mix thoroughly by pipetting:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Ultra II End-prep reaction buffer <ul style="list-style-type: none"> - Volume per reaction: 1.75 µl - For X24 samples: 52.5 µl - For X48 samples: 105 µl - For X96 samples: 210 µl <input type="checkbox"/> Ultra II End-prep enzyme mix <ul style="list-style-type: none"> - Volume per reaction: 0.75 µl - For X24 samples: 22.5 µl - For X48 samples: 45 µl - For X96 samples: 90 µl <p><input type="checkbox"/> Using a stepper pipette or multi-channel pipette, add 2.5 µl of the end-prep master mix to each well containing 12.5 µl sample.</p> <p><input type="checkbox"/> Ensure the reactions are thoroughly mixed by pipetting. Seal the End-prep plate and spin down briefly.</p> <p><input type="checkbox"/> Using a thermal cycler, incubate the plate at 20°C for 5 minutes and 65°C for 5 minutes.</p>	
<p>Take forward the end-prepped DNA into the native barcode ligation step.</p>	
<p>Native barcode ligation</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> To monitor cross-contamination events, we recommend that the negative control is carried through this process and a barcode is used to sequence this control.</p>	
<p>Thaw Native Barcodes at RT. Use one barcode per sample.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Spin down barcodes before opening tubes/piercing plates. <input type="checkbox"/> Pipette mix the entire content of each barcode 10 times before use. <input type="checkbox"/> Once thawed, keep the barcodes on ice. <p>Prepare the NEB Blunt/TA Ligase Master Mix according to the manufacturer's instructions, and place on ice:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw the reagents at RT. <input type="checkbox"/> Spin down the reagent tubes for 5 seconds. <input type="checkbox"/> Ensure the reagents are fully mixed by performing 10 full volume pipette mixes. <p><input type="checkbox"/> Thaw the Short Fragment Buffer (SFB) at RT and mix by vortexing. Then spin down and place on ice.</p> <p>In a clean 96-well plate (Native barcode ligation plate), add the following reagents in the following order per well:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 3 µl Nuclease-free water <input type="checkbox"/> 0.75 µl End-prepped DNA <input type="checkbox"/> 1.25 µl Native Barcode <input type="checkbox"/> 5 µl Blunt/TA Ligase Master Mix <p><input type="checkbox"/> Mix contents thoroughly by pipetting, seal the Native barcode ligation plate and spin down briefly.</p> <p><input type="checkbox"/> Using a thermal cycler, incubate the Native barcode ligation plate at 20°C for 20 mins and at 65°C for 10 mins.</p>	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Pool all the barcoded samples in a clean 1.5 ml Eppendorf DNA LoBind tube, checking the base of the plate to ensure all the liquid has been pooled. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Add a 0.4X volume of AMPure XP beads to the pooled reaction and mix thoroughly by pipetting: <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 80% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet the beads on a magnet for 5 minutes. Keep the tube on the magnet until the eluate is clear, and pipette off the supernatant. <input type="checkbox"/> Wash the beads by adding 700 µl Short Fragment Buffer (SFB). Flick the beads to resuspend, then return the tube to the magnetic rack and allow the beads to pellet. Keep the tube on the magnet until the eluate is clear and colourless. Remove the supernatant using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Keep the tube on the magnetic rack and wash the beads with 100 µl of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Spin down and place the tube back on the magnetic rack. Pipette off any residual ethanol. Allow the pellet to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 35 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnetic rack until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 35 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	
<p>Take forward the barcoded DNA library to the adapter ligation and clean-up step. However, at this point it is also possible to store the sample at 4°C overnight.</p>	
<p>Adapter ligation and clean-up</p>	
<p>Adapter Mix II Expansion use</p> <p>Prepare the NEBNext Quick Ligation Reaction Module according to the manufacturer's instructions, and place on ice:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw the reagents at RT. <input type="checkbox"/> Spin down the reagent tubes for 5 seconds. <input type="checkbox"/> Ensure the reagents are fully mixed by performing 10 full volume pipette mixes. <p>Note: Do NOT vortex the Quick T4 DNA Ligase.</p>	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Do not vortex the Quick T4 DNA Ligase. 	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Thaw the Elution Buffer (EB) and the Short Fragment Buffer (SFB) at RT, before mixing by vortexing. Then spin down and place on ice.</p> <p><input type="checkbox"/> Spin down the Adapter Mix II (AMII), pipette mix and place on ice.</p> <p>Perform the adapter ligation of the pooled and barcoded DNA. In a clean 1.5 ml Eppendorf LoBind tube, add the reagents in the following order:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 30 µl Pooled barcoded sample <input type="checkbox"/> 5 µl Adapter Mix II (AMII) <input type="checkbox"/> 10 µl NEBNext Quick Ligation Reaction Buffer (5X) <input type="checkbox"/> 5 µl Quick T4 DNA Ligase <p><input type="checkbox"/> Ensure the components are thoroughly mixed by pipetting, and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> The next clean-up step uses Short Fragment Buffer (SFB) and not 80% ethanol to wash the beads. The use of ethanol will be detrimental to the sequencing reaction.</p>	
<p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Add 20 µl of resuspended AMPure XP beads to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p> <p><input type="checkbox"/> Spin down the sample and pellet the beads on a magnet for 5 minutes. Keep the tube on the magnet until the eluate is clear and colourless, and pipette off the supernatant.</p> <p><input type="checkbox"/> Wash the beads by adding 125 µl Short Fragment Buffer (SFB). Flick the beads to resuspend, then return the tube to the magnetic rack and allow the beads to pellet. Keep the tube on the magnet until the eluate is clear and colourless. Remove the supernatant using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet by pipetting in 15 µl Elution Buffer (EB). Spin down and incubate for 5 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute.</p> <p><input type="checkbox"/> Remove and retain 15 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube.</p>	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> We recommend loading 50 fmol of final prepared library onto a flow cell.</p>	
<p>The prepared library is used for loading onto the flow cell. Store the library on ice until ready to load.</p>	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Priming and loading the SpotON flow cell</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT before mixing the reagents by vortexing, and spin down at RT. <input type="checkbox"/> To prepare the flow cell priming mix, add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing at RT. <input type="checkbox"/> Open the MinION device lid and slide the flow cell under the clip. <input type="checkbox"/> Slide the flow cell priming port cover clockwise to open the priming port. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores. 	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes. During this time, prepare the library for loading by following the steps below. <input type="checkbox"/> Thoroughly mix the contents of the Loading Beads (LB) by pipetting. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> The Loading Beads (LB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use. 	
<p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 37.5 µl Sequencing Buffer (SQB) <input type="checkbox"/> 25.5 µl Loading Beads (LB), mixed immediately before use <input type="checkbox"/> 12 µl DNA library <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. 	

Ligation sequencing influenza whole genome (SQK-LSK109 with EXP-NBD196)

Version: INF_9166_v109_revB_24Aug2022
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Add 75 µl of the prepared library to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION device lid.	
<p>IMPORTANT</p> <input type="checkbox"/> Required settings in MinkNOW	
<p>Flow cell reuse and returns</p>	
<input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C. <input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore.	
<p>IMPORTANT</p> <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.	