

Ligation sequencing gDNA - native barcoding (SQK-LSK109 with EXP-NBD104 and EXP-NBD114)

Version: NBE_9065_v109_revAP_14Aug2019
Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 1 µg (or 100-200 fmol) high molecular weight genomic DNA for every sample to be barcoded	<input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (NEB, M0367)	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
<input type="checkbox"/> 1.5-3 µg (or 150-300 fmol) high molecular weight genomic DNA for every sample to be barcoded, if using R10.3 flow cells	<input type="checkbox"/> NEBNext FFPE Repair Mix (NEB, M6630)	<input type="checkbox"/> Magnetic rack, suitable for 1.5 ml Eppendorf tubes
<input type="checkbox"/> OR 100+ ng high molecular weight genomic DNA if performing DNA fragmentation	<input type="checkbox"/> NEBNext Ultra II End repair/dA-tailing Module (NEB, E7546)	<input type="checkbox"/> Microfuge
<input type="checkbox"/> Native Barcoding Expansion 1-12 (EXP-NBD104) and 13-24 (EXP-NBD114) if multiplexing more than 12 samples	<input type="checkbox"/> NEBNext Quick Ligation Module (NEB, E6056)	<input type="checkbox"/> Vortex mixer
<input type="checkbox"/> Ligation Sequencing Kit (SQK-LSK109)	<input type="checkbox"/> Agencourt AMPure XP beads (Beckman Coulter, A63881)	<input type="checkbox"/> Thermal cycler
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Ice bucket with ice
<input type="checkbox"/> Adapter Mix II Expansion (EXP-AMII001)	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Timer
	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, AM9937)	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water	
INSTRUCTIONS	NOTES/OBSERVATIONS	
<p>DNA repair and end-prep</p> <p>Prepare the NEBNext FFPE DNA Repair Mix and NEBNext Ultra II End Repair / dA-tailing Module reagents in accordance with manufacturer's instructions, and place on ice.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw all reagents on ice. <input type="checkbox"/> Flick and/or invert the reagent tubes to ensure they are well mixed. Note: Do not vortex the FFPE DNA Repair Mix or Ultra II End Prep Enzyme Mix. <input type="checkbox"/> Always spin down tubes before opening for the first time each day. <input type="checkbox"/> The Ultra II End Prep Buffer and FFPE DNA Repair Buffer may have a little precipitate. Allow the mixture to come to RT and pipette the buffer up and down several times to break up the precipitate, followed by vortexing the tube for 30 seconds to solubilise any precipitate. Note: It is important the buffers are mixed well by vortexing. <input type="checkbox"/> The FFPE DNA Repair Buffer may have a yellow tinge and is fine to use if yellow. <p>Prepare the DNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <input type="checkbox"/> For R9.4.1 flow cells, transfer 1 µg (or 100-200 fmol) genomic DNA into a 1.5 ml Eppendorf DNA LoBind tube, or 1.5-3 µg (or 150-300 fmol) genomic DNA if using R10.3 flow cells. <input type="checkbox"/> Adjust the volume to 48 µl with Nuclease-free water <input type="checkbox"/> Mix thoroughly by flicking the tube <input type="checkbox"/> Spin down briefly in a microfuge 		

Ligation sequencing gDNA - native barcoding (SQK-LSK109 with EXP-NBD104 and EXP-NBD114)

Version: NBE_9065_v109_revAP_14Aug2019
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>In a 0.2 ml thin-walled PCR tube, mix the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 48 µl DNA <input type="checkbox"/> 3.5 µl NEBNext FFPE DNA Repair Buffer <input type="checkbox"/> 3.5 µl Ultra II End-prep reaction buffer <input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix <input type="checkbox"/> 2 µl NEBNext FFPE DNA Repair Mix <p><input type="checkbox"/> Mix well by pipetting using wide-bore pipette tips. Alternatively, if you are concerned about preserving the integrity of very long DNA fragments, mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Using a thermal cycler, incubate at 20°C for 5 minutes and 65°C for 5 minutes.</p> <p><input type="checkbox"/> Transfer the DNA sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by flicking the tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet until supernatant is clear and colourless. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 25 µl Nuclease-free water. Spin down and incubate for 2 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute.</p> <p><input type="checkbox"/> Remove and retain 25 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer - recovery aim >700 ng.</p>	
<p>Take forward the repaired and end-prepped DNA into the native barcode ligation step. However, at this point it is also possible to store the sample at 4°C overnight.</p>	
<p>Native barcode ligation</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the native barcodes at RT. Use one barcode per sample. Individually mix the barcodes by pipetting, spin down, and place them on ice. <input type="checkbox"/> Select a unique barcode for every sample to be run together on the same flow cell, from the provided 24 barcodes. Up to 24 samples can be barcoded and combined in one experiment. 	

Ligation sequencing gDNA - native barcoding (SQK-LSK109 with EXP-NBD104 and EXP-NBD114)

Version: NBE_9065_v109_revAP_14Aug2019
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Dilute 500 ng (750 ng if sequencing on R10.3 flow cells) of each end-prepped sample to be barcoded to 22.5 µl in Nuclease-free water.</p> <p>Add the reagents in the order given below, mixing by flicking the tube between each sequential addition:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 22.5 µl 500 ng end-prepped DNA (750 ng if using R10.3 flow cells) <input type="checkbox"/> 2.5 µl Native Barcode <input type="checkbox"/> 25 µl Blunt/TA Ligase Master Mix <p><input type="checkbox"/> Mix well by pipetting using wide-bore pipette tips. Alternatively, if you are concerned about preserving the integrity of very long DNA fragments, mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</p> <p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Add 50 µl of resuspended AMPure XP beads to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant when clear and colourless.</p> <p><input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 26 µl Nuclease-free water. Incubate for 2 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 26 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube.</p>	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	
<p><input type="checkbox"/> Pool equimolar amounts of each barcoded sample into a 1.5 ml Eppendorf DNA LoBind tube, ensuring that sufficient sample is combined to produce a pooled sample of 700 ng total (1050 ng if sequencing on R10.3 flow cells).</p> <p><input type="checkbox"/> Quantify 1 µl of pooled and barcoded DNA using a Qubit fluorometer.</p> <p><input type="checkbox"/> Dilute 700 ng (1050 ng for R10.3 flow cells) pooled sample to 65 µl in Nuclease-free water.</p>	
<p>Take forward the pooled samples into the next step. However, at this point it is also possible to store the sample at 4°C overnight.</p>	

Ligation sequencing gDNA - native barcoding (SQK-LSK109 with EXP-NBD104 and EXP-NBD114)

Version: NBE_9065_v109_revAP_14Aug2019
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Adapter ligation and clean-up</p> <p>Adapter Mix II Expansion use</p> <p><input type="checkbox"/> Thaw the Elution Buffer (EB) and NEBNext Quick Ligation Reaction Buffer (5x) at RT, mix by vortexing, spin down and place on ice. Check the contents of each tube are clear of any precipitate.</p> <p><input type="checkbox"/> Spin down the T4 Ligase and the Adapter Mix II (AMII), and place on ice.</p> <p>IMPORTANT</p> <p>Depending on the wash buffer (LFB or SFB) used, the clean-up step after adapter ligation is designed to either enrich for DNA fragments of >3 kb, or purify all fragments equally.</p> <p><input type="checkbox"/> To enrich for DNA fragments of 3 kb or longer, use Long Fragment Buffer (LFB)</p> <p><input type="checkbox"/> To retain DNA fragments of all sizes, use Short Fragment Buffer (SFB)</p> <p><input type="checkbox"/> To enrich for DNA fragments of 3 kb or longer, thaw one tube of Long Fragment Buffer (LFB) at RT, mix by vortexing, spin down and place on ice.</p> <p><input type="checkbox"/> To retain DNA fragments of all sizes, thaw one tube of Short Fragment Buffer (SFB) at RT, mix by vortexing, spin down and place on ice.</p> <p>Taking the pooled and barcoded DNA, perform adapter ligation as follows, mixing by flicking the tube between each sequential addition.</p> <ul style="list-style-type: none"> <input type="checkbox"/> 65 µl 700 ng (1050 ng if using R10.3 flow cells) pooled barcoded sample <input type="checkbox"/> 5 µl Adapter Mix II (AMII) <input type="checkbox"/> 20 µl NEBNext Quick Ligation Reaction Buffer (5X) <input type="checkbox"/> 10 µl Quick T4 DNA Ligase <p><input type="checkbox"/> Ensure the components are thoroughly mixed by pipetting, and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</p> <p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Add 50 µl of resuspended AMPure XP beads to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Place on a magnetic rack, allow beads to pellet and pipette off supernatant.</p> <p><input type="checkbox"/> Wash the beads by adding either 250 µl Long Fragment Buffer (LFB) or 250 µl Short Fragment Buffer (SFB). Flick the beads to resuspend, spin down, then return the tube to the magnetic rack and allow the beads to pellet. Remove the supernatant using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 15 µl Elution Buffer (EB). Spin down and incubate for 10 minutes at RT. For high molecular weight DNA, incubating at 37°C can improve the recovery of long fragments.</p>	

Ligation sequencing gDNA - native barcoding (SQK-LSK109 with EXP-NBD104 and EXP-NBD114)

Version: NBE_9065_v109_revAP_14Aug2019
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute. <input type="checkbox"/> Remove and retain 15 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Quantify 1 µl of adapter ligated and barcoded DNA using a Qubit fluorometer - recovery aim ~430 ng. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> We recommend loading 5-50 fmol of final prepared library onto a flow cell. 	
<p>The prepared library is used for loading onto the flow cell. Store the library on ice until ready to load.</p>	
<p></p>	
<p>Priming and loading the SpotON flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT before mixing the reagents by vortexing, and spin down at RT. <input type="checkbox"/> To prepare the flow cell priming mix, add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing at RT. <input type="checkbox"/> Open the MinION device lid and slide the flow cell under the clip. <input type="checkbox"/> Slide the flow cell priming port cover clockwise to open the priming port. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores. 	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes. During this time, prepare the library for loading by following the steps below. <input type="checkbox"/> Thoroughly mix the contents of the Loading Beads (LB) by pipetting. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> The Loading Beads (LB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use. 	

Ligation sequencing gDNA - native barcoding (SQK-LSK109 with EXP-NBD104 and EXP-NBD114)

Version: NBE_9065_v109_revAP_14Aug2019
 Last update: 10/03/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 37.5 µl Sequencing Buffer (SQB) <input type="checkbox"/> 25.5 µl Loading Beads (LB), mixed immediately before use <input type="checkbox"/> 12 µl DNA library <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of the prepared library to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION device lid. 	
<p>Flow cell reuse and returns</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C. <input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol. 	