

# 1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> <1 µg of each DNA sample to be barcoded in 45 µl	<input type="checkbox"/> Agencourt AMPure XP beads	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
<input type="checkbox"/> PCR Barcoding Expansion Pack 1-96 (EXP-PBC096)	<input type="checkbox"/> NEBNext® Companion Module for Oxford Nanopore Technologies® Ligation Sequencing (cat # E7180S). Alternatively, you can use the three NEBNext® products below:	<input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes
<input type="checkbox"/> Ligation Sequencing Kit (SQK-LSK109)	<input type="checkbox"/> NEBNext FFPE Repair Mix (M6630)	<input type="checkbox"/> Microfuge
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> NEBNext End repair / dA-tailing Module (E7546)	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> NEBNext Quick Ligation Module (E6056)	<input type="checkbox"/> Thermal cycler
	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Ice bucket with ice
	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Timer
	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937)	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water	
	<input type="checkbox"/> LongAmp Taq 2X Master Mix (e.g. NEB M0287)	

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Preparing input DNA</b></p> <p><input type="checkbox"/> Record the quality, quantity and size of the DNA.</p> <p><b>IMPORTANT</b></p> <p>Criteria for input DNA</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Purity as measured using Nanodrop - OD 260/280 of 1.8 and OD 260/230 of 2.0-2.2</li> <li><input type="checkbox"/> Input mass, as measured by Qubit - &lt;1 µg, or 100-200 fmol for short-fragment libraries</li> <li><input type="checkbox"/> No detergents or surfactants in the buffer</li> </ul> <p>Prepare the DNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Transfer &lt;1 µg DNA into a DNA LoBind tube</li> <li><input type="checkbox"/> Adjust the volume to 45 µl with Nuclease-free water</li> <li><input type="checkbox"/> Mix thoroughly by inversion avoiding unwanted shearing</li> <li><input type="checkbox"/> Spin down briefly in a microfuge</li> </ul>	

1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Check your flow cell</b></p> <p><input type="checkbox"/> Set up the MinION, flow cell and host computer</p> <p>Once successfully plugged in, you will see a light and hear the fan.</p> <p>Open the MinkNOW GUI from the desktop icon and establish a local or remote connection.</p> <p><input type="checkbox"/> Plug the MinION into the computer.</p> <p><input type="checkbox"/> Choose the flow cell type from the selector box. Then check the "Available" box.                      Note: if you are using flow cells from your Starter Pack, please select FLO-MIN106.</p> <p>Click "Check flow cells" at the bottom of the screen.</p> <p><input type="checkbox"/> R9.4.1 FLO-MIN106</p> <p><input type="checkbox"/> R9.5.1 FLO-MIN107</p> <p><input type="checkbox"/> Click "Start test".</p> <p><input type="checkbox"/> Check the number of active pores available for the experiment, reported in the System History panel when the check is complete.</p>	
<p>Flow cell check complete.</p>	
<p><b>End-prep</b></p> <p>In a 0.2 ml 96 well PCR plate, set up the end-repair / dA-tailing reactions as follows:</p> <p><input type="checkbox"/> 45 µl &lt;1 µg DNA</p> <p><input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer</p> <p><input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix</p> <p><input type="checkbox"/> 5 µl Nuclease-free water</p> <p><input type="checkbox"/> Mix by pipetting.</p> <p><input type="checkbox"/> Seal the plate with adhesive film or PCR strip caps, spin down in a centrifuge and incubate for 5 minutes at 20 °C and 5 minutes at 65 °C using the thermal cycler.</p> <p><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</p> <p><input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting.</p> <p><input type="checkbox"/> Allow DNA to bind to beads for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare sufficient fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Place on magnetic rack, allow beads to pellet and pipette off supernatant.</p> <p><input type="checkbox"/> Keep on magnet, wash beads with 180 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p>	

1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <li><input type="checkbox"/> Cover the plate with adhesive film and leave plate on magnet for 2 minutes to allow residual liquid to collect at the bottom. Remove the adhesive film, return the plate to the magnet and aspirate residual wash solution.</li> <li><input type="checkbox"/> Briefly incubate the plate on a thermal cycler at 37° C with the lid open and the plate wells unsealed.</li> <li><input type="checkbox"/> Remove the plate from the magnet and resuspend pellet in 31 µl Nuclease-free water. Incubate for 2 minutes at RT.</li> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</li> <li><input type="checkbox"/> Remove eluate once it is clear and colourless. Transfer each eluted sample to a new 96-well PCR plate.</li> <li><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer - recovery aim &gt; 700 ng.</li> </ul>	
<p>Take forward approximately 700 ng of end-prepped DNA in 30 µl into adapter ligation.</p>	
<p><b>Ligation of Barcode Adapter</b></p>	
<p>Add the reagents to a fresh 96 well plate, in the order given below:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 30 µl End prep DNA</li> <li><input type="checkbox"/> 20 µl Barcode Adapter</li> <li><input type="checkbox"/> 50 µl Blunt/TA Ligase Master Mix</li> <li><input type="checkbox"/> Mix by pipetting.</li> <li><input type="checkbox"/> Seal the plate with adhesive film or PCR strip caps and briefly spin down in a plate spinner.</li> <li><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</li> <li><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</li> <li><input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube.</li> <li><input type="checkbox"/> Allow DNA to bind to beads for 5 minutes at RT.</li> <li><input type="checkbox"/> Prepare sufficient fresh 70% ethanol in Nuclease-free water.</li> <li><input type="checkbox"/> Place on magnetic rack, allow beads to pellet and pipette off supernatant.</li> <li><input type="checkbox"/> Keep on magnet, wash beads with 180 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</li> <li><input type="checkbox"/> Repeat the previous step.</li> <li><input type="checkbox"/> Cover the plate with adhesive film and leave plate on magnet for 2 minutes to allow residual liquid to collect at the bottom. Remove the adhesive film, return the plate to the magnet and aspirate residual wash solution.</li> <li><input type="checkbox"/> Briefly incubate the plate on a thermal cycler at 37° C with the lid open and the plate wells unsealed.</li> <li><input type="checkbox"/> Remove the plate from the magnet and resuspend pellet in 25 µl Nuclease-free water. Incubate for 2 minutes at RT.</li> </ul>	


1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</li> <li><input type="checkbox"/> Remove eluate once it is clear and colourless. Transfer each eluted sample to a new 96-well PCR plate.</li> <li><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer.</li> <li><input type="checkbox"/> Dilute the library to a concentration of 10 ng/µl with Nuclease-free water or 10 mM Tris-HCl pH 8.5</li> </ul>	
<p><b>Barcoding PCR</b></p>	
<p>Set up a barcoding PCR reaction as follows for each library:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 1 µl PCR Barcode (one of BC1-BC96)</li> <li><input type="checkbox"/> 2 µl 10 ng/µl adapter ligated template</li> <li><input type="checkbox"/> 25 µl LongAmp Taq 2x master mix</li> <li><input type="checkbox"/> 22 µl Nuclease-free water</li> </ul> <p><input type="checkbox"/> Mix by pipetting.</p> <p><input type="checkbox"/> Seal the plate with adhesive film or PCR strip caps and briefly spin down in a plate spinner.</p> <p>Amplify using the following cycling conditions:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Initial denaturation 3 mins @ 95 °C (1 cycle)</li> <li><input type="checkbox"/> Denaturation 15 secs @ 95 °C (15-18 (b) cycles)</li> <li><input type="checkbox"/> Annealing 15 secs (a) @ 62 °C (a) (15-18 (b) cycles)</li> <li><input type="checkbox"/> Extension dependent on length of target fragment (d) @ 65 °C (c) (15-18 (b) cycles)</li> <li><input type="checkbox"/> Final extension dependent on length of target fragment (d) @ 65 °C (1 cycle)</li> <li><input type="checkbox"/> Hold @ 4 °C</li> </ul> <p><input type="checkbox"/> Purify the barcoded DNA using standard methods which are suitable for the fragment size.</p> <p><input type="checkbox"/> Quantify the barcoded library using standard techniques, and pool all barcoded libraries in the desired ratios in a 1.5 ml DNA LoBind Eppendorf tube.</p> <p><input type="checkbox"/> Prepare 1 µg of pooled barcoded libraries in 47 µl Nuclease-free water.</p>	
<p>This pooled library is now ready to be end-repaired and adapted for nanopore sequencing.</p>	
<p><b>DNA repair and end-prep</b></p>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Thaw DNA CS (DCS) at RT, spin down, mix by pipetting, and place on ice.</li> <li><input type="checkbox"/> Prepare the NEBNext FFPE DNA Repair Mix and NEBNext End repair / dA-tailing Module reagents in accordance with manufacturer's instructions, and place on ice.</li> </ul> <p>In a 0.2 ml thin-walled PCR tube, mix the following:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 1 µl DNA CS</li> <li><input type="checkbox"/> 47 µl DNA</li> <li><input type="checkbox"/> 3.5 µl NEBNext FFPE DNA Repair Buffer</li> <li><input type="checkbox"/> 2 µl NEBNext FFPE DNA Repair Mix</li> <li><input type="checkbox"/> 3.5 µl Ultra II End-prep reaction buffer</li> <li><input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix</li> </ul>	

1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <li><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</li> <li><input type="checkbox"/> Using a thermal cycler, incubate at 20° C for 5 minutes and 65° C for 5 mins.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> AMPure XP bead clean-up</li> </ul>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</li> <li><input type="checkbox"/> Transfer the DNA sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</li> <li><input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by flicking the tube.</li> <li><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</li> <li><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</li> <li><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</li> <li><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</li> <li><input type="checkbox"/> Repeat the previous step.</li> <li><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</li> <li><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 61 µl Nuclease-free water. Incubate for 2 minutes at RT.</li> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</li> <li><input type="checkbox"/> Remove and retain 61 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</li> </ul>	
<p>Take forward the repaired and end-prepped DNA into the adapter ligation step. However, at this point it is also possible to store the sample at 4° C overnight.</p>	
<p><b>Adapter ligation and clean-up</b></p>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Although the recommended 3rd party ligase is supplied with its own buffer, the ligation efficiency of Adapter Mix (AMX) is higher when using Ligation Buffer supplied within SQK-LSK109.</li> </ul>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Spin down Adapter Mix (AMX) and T4 Ligase from the NEBNext Quick Ligation Module (E6056), and place on ice.</li> <li><input type="checkbox"/> Thaw Ligation Buffer (LNB) at RT, spin down and mix by pipetting. Due to viscosity, vortexing this buffer is ineffective. Place on ice immediately after thawing and mixing.</li> <li><input type="checkbox"/> Thaw the Elution Buffer (EB) at RT, mix by vortexing, spin down and place on ice.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Depending on the wash buffer used in this section, the clean-up step after adapter ligation is designed to either enrich for long DNA fragments, or purify all fragments equally.</li> </ul>	

1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> To enrich for DNA fragments of 3 kb or longer, thaw one tube of L Fragment Buffer (LFB) at RT, mix by vortexing, spin down and place on ice.</p> <p><input type="checkbox"/> To retain DNA fragments shorter than 3 kb (by purifying fragments of all sizes), thaw one tube of S Fragment Buffer (SFB) at RT, mix by vortexing, spin down and place on ice.</p> <p>In a 1.5 ml Eppendorf DNA LoBind tube, mix in the following order:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 60 µl DNA sample from the previous step</li> <li><input type="checkbox"/> 25 µl Ligation Buffer (LNB)</li> <li><input type="checkbox"/> 10 µl NEBNext Quick T4 DNA Ligase</li> <li><input type="checkbox"/> 5 µl Adapter Mix (AMX)</li> </ul> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> If you have omitted the AMPure purification step after DNA repair and end-prep, incubating the ligation reaction for longer than 10 minutes is not recommended.</p>	
<p><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</p> <p><input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Wash the beads by adding either 250 µl Long Fragment Buffer (LFB) or 250 µl S Fragment Buffer (SFB). Flick the beads to resuspend, then return the tube to the magnetic rack and allow the beads to pellet. Remove the supernatant using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 15 µl Elution Buffer (EB). Incubate for 10 minutes at RT. For high molecular weight DNA, incubating at 37° C can improve the recovery of long fragments.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p>Remove and retain 15 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Remove and retain the eluate which contains the DNA library in a clean 1.5 ml Eppendorf DNA LoBind tube</li> <li><input type="checkbox"/> Dispose of the pelleted beads</li> </ul>	
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	
<p>The prepared library is used for loading into the flow cell. Store the library on ice until ready to load.</p>	

1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Recommendation on library loading</p>	
<p><b>Priming and loading the SpotON flow cell</b></p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Please note that the Sequencing Tether (SQT) tube will NOT be used in this protocol. It is provided in the kit for potential future product compatibility.</p>	
<p><input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FLB) at RT before placing the tubes on ice as soon as thawing is complete.</p> <p><input type="checkbox"/> Mix the Sequencing Buffer (SQB) and Flush Buffer (FLB) tubes by vortexing, spin down and return to ice.</p> <p><input type="checkbox"/> Spin down the Flush Tether (FLT) tube, mix by pipetting, and return to ice.</p> <p><input type="checkbox"/> Open the lid of the nanopore sequencing device and slide the flow cell's priming port cover clockwise so that the priming port is visible.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Care must be taken when drawing back buffer from the flow cell. The array of pores must be covered by buffer at all times. Removing more than 20-30 µl risks damaging the pores in the array.</p>	
<p>After opening the priming port, check for small bubble under the cover. Draw back a small volume to remove any bubble (a few µl):</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 µl</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p><input type="checkbox"/> Prepare the flow cell priming mix: add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by pipetting up and down.</p> <p><input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes.</p> <p><input type="checkbox"/> Thoroughly mix the contents of the Loading Beads (LB) by pipetting.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> The Loading Beads (LB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use.</p>	
<p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 37.5 µl Sequencing Buffer (SQB)</li> <li><input type="checkbox"/> 25.5 µl Loading Beads (LB), mixed immediately before use</li> <li><input type="checkbox"/> 12 µl DNA library</li> </ul> <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li> <li><input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li> </ul>	

Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <li><input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading.</li> <li><input type="checkbox"/> Add 75 µl of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</li> <li><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION lid.</li> </ul>	
<p><b>Starting a sequencing run</b></p>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Double-click the MinKNOW icon located on the desktop to open the MinKNOW GUI.</li> <li><input type="checkbox"/> If your MinION was disconnected from the computer, plug it back in.</li> <li><input type="checkbox"/> Choose the flow cell type from the selector box. Then check the "Available" box.</li> <li><input type="checkbox"/> Click the "New Experiment" button at the bottom left of the GUI.</li> </ul> <p>On the New experiment popup screen, select the running parameters for your experiment from the individual tabs.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Fast - This model is able to keep up with a full experiment on a GridION or MinIT. It implements a lite version of the Flip-flop algorithm. This delivers the same level of accuracy as the previous default transducer models, which have been deployed prior to the introduction of Flip-flop.</li> <li><input type="checkbox"/> HAC - The High accuracy (HAC) Flip-flop model provides a higher consensus/raw read accuracy than the 'Fast' model. It contains a more computationally intense Flip-flop architecture that can deliver higher accuracy using the same data produced by nanopore sequencing. It is currently 5-8 times slower than the fast model so users should ensure their data transfer, disk-space and device utilisation is scaled appropriately for this.</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Click "Start run".</li> </ul> <p>Allow the script to run to completion.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> The MinKNOW Experiment page will indicate the progression of the script; this can be accessed through the "Experiment" tab that will appear at the top right of the screen</li> <li><input type="checkbox"/> Monitor messages in the Message panel in the MinKNOW GUI</li> </ul> <p>Basecalled read files</p>	
<p><b>Progression of MinKNOW protocol script</b></p>	
<p>The running experiment screen</p> <p>Experiment summary information</p> <p>Check the number of active pores reported in the MUX scan are similar (within 10-15%) to those reported at the end of the Flow Cell Check</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> If there is a significant reduction in the numbers, restart MinKNOW.</li> <li><input type="checkbox"/> If the numbers are still significantly different, close down the host computer and reboot.</li> <li><input type="checkbox"/> When the numbers are similar to those reported at the end of the Flow Cell Check, restart the experiment on the Connection page. There is no need to load any additional library after restart.</li> <li><input type="checkbox"/> Stopping the experiment is achieved by clicking "Stop run" button at the top of the screen.</li> </ul> <p>Data acquisition will stop, but the software will continue basecalling unless the user clicks the "Stop basecalling" button.</p>	



1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
 Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Check that the temperature has reached 34° C.</p> <p>Check pore occupancy in the channel panel at the top of the experimental view.</p> <p><input type="checkbox"/> A good library will be indicated by a higher proportion of light green channels in Sequencing than are in Pore. The combination of Sequencing and Pore indicates the number of active pores at any point in time. A low proportion of Sequencing channels will reduce the throughput of the run.</p> <p><input type="checkbox"/> Recovering indicates channels that may become available for sequencing again. A high proportion of this may indicate additional clean up steps are required during your library preparation.</p> <p><input type="checkbox"/> Inactive indicates channels that are no longer available for sequencing. A high proportion of these as soon as the run begins may indicate an osmotic imbalance.</p> <p><input type="checkbox"/> Unclassified are channels that have not yet been assigned one of the above classifications</p> <p><input type="checkbox"/> Monitor the pore occupancy</p> <p><input type="checkbox"/> Monitor the disk space</p> <p>Duty time plots</p> <p><input type="checkbox"/> Monitor the development of the read length histogram.</p> <p>Cumulative output</p> <p><input type="checkbox"/> the number of reads that have been sequenced and basecalled; and whether the reads have passed of failed the quality filters</p> <p><input type="checkbox"/></p> <p>MinKNOW 19.05 has several new interpretations and graphics of the data obtained from a sequencing run.</p> <p>Temperature and Bias voltage graph</p> <p>Translocation speed and Qscore graphs</p> <p>Cumulative output of multiple flow cells - experiment view on PromethION</p> <p>Trace viewer</p>	
<p><b>Onward analysis of MinKNOW basecalled data</b></p>	
<p><input type="checkbox"/> Open the Desktop Agent using the desktop shortcut.</p> <p><input type="checkbox"/> Click on the New Workflow tab in the Desktop Agent and select the FASTQ barcoding workflow.</p> <p>Select the workflow parameters.</p> <p><input type="checkbox"/> Select the quality score cut-off</p> <p><input type="checkbox"/> Check the correct settings are selected in the Desktop Agent.</p> <p><input type="checkbox"/> Click "Start Run" to start data analysis.</p> <p><input type="checkbox"/> Follow the progression of upload and download of read files in the Desktop Agent.</p>	

1D PCR barcoding (96) genomic DNA (SQK-LSK109)

Version: PBGE96\_9068\_v109\_revF\_23May2018  
Last update: 02/05/2019



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Click on VIEW REPORT.</p> <ul style="list-style-type: none"><li><input type="checkbox"/> Click on VIEW REPORT to navigate to the EPI2ME website; this can be done at any point during data exchange</li><li><input type="checkbox"/> Return to the Desktop Agent to see progression of the exchange</li></ul>	
<p><b>Close down MinKNOW and the Desktop Agent</b></p>	
<ul style="list-style-type: none"><li><input type="checkbox"/> Quit Desktop Agent using the close x.</li><li><input type="checkbox"/> Quit MinKNOW by closing down the web GUI.</li><li><input type="checkbox"/> Disconnect the MinION.</li></ul>	
<p><b>Prepare the flow cell for re-use or return to Oxford Nanopore.</b></p>	
<ul style="list-style-type: none"><li><input type="checkbox"/> If you would like to reuse the flow cell, follow the Wash Kit instructions and store the washed flow cell at 2-8 °C, OR</li><li><input type="checkbox"/> Follow the returns procedure by washing out the MinION Flow Cell ready to send back to Oxford Nanopore.</li></ul>	