

# 1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> PCR Barcoding Expansion 1-12 (EXP-PBC001)	<input type="checkbox"/> NEBNext End repair / dA-tailing Module (E7546)	<input type="checkbox"/> Thermal cycler
<input type="checkbox"/> Ligation Sequencing Kit 1D (SQK-LSK108)	<input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (M0367)	<input type="checkbox"/> Microfuge
<input type="checkbox"/> Library Loading Bead Kit (EXP-LLB001)	<input type="checkbox"/> Agencourt AMPure XP beads	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water	<input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes
	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Heating block at 37° C capable of taking 1.5 ml tubes
	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Pipettes P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937)	
	<input type="checkbox"/> LongAmp Taq 2X Master Mix (e.g. NEB M0287)	
	<input type="checkbox"/> Pipette tips P2, P10, P20, P100, P200, P1000	
INSTRUCTIONS		NOTES/OBSERVATIONS
<b>Preparing input DNA</b>		
Prepare the DNA in Nuclease-free water. <ul style="list-style-type: none"> <li><input type="checkbox"/> Transfer 1-1.5 µg genomic DNA into a DNA LoBind tube</li> <li><input type="checkbox"/> Adjust the volume to 46 µl with Nuclease-free water</li> <li><input type="checkbox"/> Mix thoroughly by inversion avoiding unwanted shearing</li> <li><input type="checkbox"/> Spin down briefly in a microfuge</li> </ul> <input type="checkbox"/> Record the quality, quantity and size of the DNA.		
<b>IMPORTANT</b>		
Criteria for input DNA <ul style="list-style-type: none"> <li><input type="checkbox"/> Purity as measured using Nanodrop - OD 260/280 of 1.8 and OD 260/230 of 2.0-2.2</li> <li><input type="checkbox"/> Input mass, as measured by Qubit - &lt;1 µg, or 100-200 fmol for short fragment libraries</li> <li><input type="checkbox"/> No detergents or surfactants in the buffer</li> </ul>		
<b>Check your flow cell</b>		
<input type="checkbox"/> Set up the MinION, flow cell and host computer		
Once successfully plugged in, you will see a light and hear the fan.		

# 1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Open the MinKNOW GUI from the desktop icon and establish a local or remote connection.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> If running a MinION on the same host computer, plug the MinION into the computer.</li> <li><input type="checkbox"/> If running a MinION on a remote computer, first enter the name or IP address of the remote host under Connect to a remote computer (if running from the Connection page), or Connections (if running from the homepage) and click Connect.</li> <li><input type="checkbox"/> Choose the flow cell type from the selector box. Then mark the flow cell as "Selected":</li> </ul> <p>Click "Check flow cells" at the bottom of the screen.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> R9.4.1 FLO-MIN106</li> <li><input type="checkbox"/> R9.5.1 FLO-MIN107</li> <li><input type="checkbox"/> Click "Start test".</li> <li><input type="checkbox"/> Check the number of active pores available for the experiment, reported in the System History panel when the check is complete.</li> </ul>	
<p>Flow cell check complete.</p>	
<p><b>DNA fragmentation</b></p>	
<p><b>OPTIONAL</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Transfer each sample of &lt;1 µg genomic DNA in 46 µl to a Covaris g-TUBE.</li> </ul> <p>Spin the g-TUBE for 1 minute at RT at the speed for the fragment size required.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Spin the g-TUBE for 1 minute</li> <li><input type="checkbox"/> Remove and check all the DNA has passed through the g-TUBE</li> <li><input type="checkbox"/> If DNA remains in the upper chamber, spin again for 1 minute at the same speed</li> </ul> <p>Invert the g-TUBE and spin again for 1 minute to collect the fragmented DNA.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Remove g-TUBE, invert the tube and replace into the centrifuge</li> <li><input type="checkbox"/> Spin the g-TUBE for 1 minute</li> <li><input type="checkbox"/> Remove and check the DNA has passed into the lower chamber</li> <li><input type="checkbox"/> If DNA remains in the upper chamber, spin again for 1 minute</li> <li><input type="checkbox"/> Remove g-TUBE</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Transfer the 46 µl fragmented DNA to a clean 1.5 ml Eppendorf DNA LoBind tube.</li> </ul>	
<p>Analyse 1 µl of the fragmented DNA for fragment size, quantity and quality.</p>	
<p>&lt;1 µg fragmented DNA in 45 µl from each sample is taken into the next step.</p>	
<p><b>End-prep</b></p>	
<p>Perform end-repair / dA-tailing of fragmented DNA as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 45 µl &lt;1 µg end-repaired DNA</li> <li><input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer</li> <li><input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix</li> <li><input type="checkbox"/> 5 µl Nuclease-free water</li> </ul>	

# 1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <li><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</li> <li><input type="checkbox"/> Transfer the sample to a 0.2 ml PCR tube, and incubate for 5 minutes at 20° C and 5 minutes at 65° C using the thermal cycler.</li> <li><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</li> <li><input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube.</li> <li><input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting.</li> <li><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</li> <li><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</li> <li><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</li> <li><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</li> <li><input type="checkbox"/> Repeat the previous step.</li> <li><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual 70% ethanol. Briefly allow to dry.</li> <li><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 31 µl Nuclease-free water. Incubate for 2 minutes at RT.</li> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</li> <li><input type="checkbox"/> Remove and retain 31 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</li> <li><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer - recovery aim &gt; 700 ng.</li> </ul>	
<p>Take forward approximately 700 ng of end-prepped DNA in 30 µl into adapter ligation.</p>	
<p><b>Ligation of Barcode Adapter</b></p>	
<p>Add the reagents in the order given below, mixing by flicking the tube between each sequential addition:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 30 µl End prep DNA</li> <li><input type="checkbox"/> 20 µl Barcode Adapter</li> <li><input type="checkbox"/> 50 µl Blunt/TA Ligase Master Mix</li> <li><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</li> <li><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</li> <li><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</li> <li><input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the reaction and mix by pipetting.</li> <li><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</li> </ul>	

# 1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual 70% ethanol. Briefly allow to dry.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 25 µl Nuclease-free water. Incubate for 2 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p>Remove and retain 15 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Remove and retain the eluate which contains the DNA library in a clean 1.5 ml Eppendorf DNA LoBind tube</p> <p><input type="checkbox"/> Dispose of the pelleted beads</p> <p><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer.</p> <p><input type="checkbox"/> Dilute the library to a concentration of 10 ng/µl with Nuclease-free water or 10 mM Tris-HCl pH 8.5</p>	
<p><b>Barcoding PCR</b></p>	
<p>Set up a barcoding PCR reaction as follows for each library:</p> <p><input type="checkbox"/> 2 µl PCR Barcode (one of BC01-BC12)</p> <p><input type="checkbox"/> 2 µl 10 ng/µl adapter ligated template</p> <p><input type="checkbox"/> 50 µl LongAmp Taq 2x master mix</p> <p><input type="checkbox"/> 46 µl Nuclease-free water</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Amplify using the following cycling conditions:</p> <p><input type="checkbox"/> Initial denaturation 3 mins @ 95 °C (1 cycle)</p> <p><input type="checkbox"/> Denaturation 15 secs @ 95 °C (12-15 (b) cycles)</p> <p><input type="checkbox"/> Annealing 15 secs (a) @ 62 °C (a) (12-15 (b) cycles)</p> <p><input type="checkbox"/> Extension dependent on length of target fragment (d) @ 65 °C (c) (12-15 (b) cycles)</p> <p><input type="checkbox"/> Final extension dependent on length of target fragment (d) @ 65 °C (1 cycle)</p> <p><input type="checkbox"/> Hold @ 4 °C</p> <p><input type="checkbox"/> Purify the barcoded DNA using standard methods which are suitable for the fragment size.</p> <p><input type="checkbox"/> Quantify the barcoded library using standard techniques, and pool all barcoded libraries in the desired ratios in a 1.5 ml DNA LoBind Eppendorf tube.</p> <p><input type="checkbox"/> Prepare 1 µg of pooled barcoded libraries in 45 µl Nuclease-free water.</p>	
<p>This pooled library is now ready to be adapted for nanopore sequencing.</p>	



# 1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>End-prep</b></p> <p>Perform end repair and dA-tailing of fragmented DNA as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 45 µl ~1 µg DNA (fragmented genomic DNA, amplicon or cDNA)</li> <li><input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer</li> <li><input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix</li> <li><input type="checkbox"/> 5 µl DNA CS</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</li> <li><input type="checkbox"/> Transfer the sample to a 0.2 ml PCR tube, and incubate for 5 minutes at 20 °C and 5 minutes at 65 °C using the thermal cycler.</li> <li><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</li> <li><input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube.</li> <li><input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting.</li> <li><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</li> <li><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</li> <li><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</li> <li><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</li> <li><input type="checkbox"/> Repeat the previous step.</li> <li><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual 70% ethanol. Briefly allow to dry.</li> <li><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 31 µl Nuclease-free water. Incubate for 2 minutes at RT.</li> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</li> <li><input type="checkbox"/> Remove and retain 31 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</li> <li><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer - recovery aim &gt; 700 ng.</li> </ul>	
<p>Take forward approximately 700 ng of end-prepped DNA in 30 µl into adapter ligation.</p>	
<p><b>Adapter ligation</b></p> <p>Thaw and prepare the kit reagents as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> ABB Buffer (ABB) at RT</li> <li><input type="checkbox"/> Elution Buffer (ELB) at RT</li> <li><input type="checkbox"/> Adapter Mix 1D (AMX1D) on ice</li> <li><input type="checkbox"/> Running Buffer with Fuel Mix (RBF) on ice</li> <li><input type="checkbox"/> Blunt/TA Ligation Master Mix on ice</li> </ul>	

1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS						
<p>Check the contents of each tube are clear of any precipitate and are thoroughly mixed before setting up the reaction.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Mix the contents of each tube by flicking</li> <li><input type="checkbox"/> Check that there is no precipitate present (DTT in the Blunt/TA Master Mix can sometimes form a precipitate)</li> <li><input type="checkbox"/> Spin down briefly before accurately pipetting the contents in the reaction</li> </ul> <p>Taking the end-prepped DNA, perform adapter ligation as follows, mixing by flicking the tube between each sequential addition.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 30 µl End-prepped DNA</li> <li><input type="checkbox"/> 20 µl Adapter Mix</li> <li><input type="checkbox"/> 50 µl Blunt/TA Ligation Master Mix</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</li> <li><input type="checkbox"/> Incubate the reaction for 10 minutes at RT.</li> </ul>							
<p><b>AMPure XP bead binding</b></p>							
<ul style="list-style-type: none"> <li><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</li> <li><input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the adapter ligation reaction from the previous step and mix by pipetting.</li> <li><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</li> <li><input type="checkbox"/> Place on magnetic rack, allow beads to pellet and pipette off supernatant.</li> <li><input type="checkbox"/> Add 140 µl of ABB Buffer (ABB) to the beads. Close the tube lid, and resuspend the beads by flicking the tube. Return the tube to the magnetic rack, allow beads to pellet and pipette off the supernatant.</li> <li><input type="checkbox"/> Repeat the previous step.</li> <li><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 15 µl Elution Buffer. Incubate for 10 minutes at RT.</li> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</li> </ul> <p>Remove and retain 15 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Remove and retain the eluate which contains the DNA library in a clean 1.5 ml Eppendorf DNA LoBind tube</li> <li><input type="checkbox"/> Dispose of the pelleted beads</li> </ul>							
<p>The prepared library is used for loading into the flow cell. Store the library on ice until ready to load.</p>							
<p><b>Before sequencing checklist</b></p> <table border="0" style="width: 100%;"> <tr> <td style="width: 33%;"><input type="checkbox"/> Prepared library on ice</td> <td style="width: 33%;"><input type="checkbox"/> Computer set up to run MinKNOW</td> <td style="width: 33%;"><input type="checkbox"/> Hardware check complete</td> </tr> <tr> <td><input type="checkbox"/> Sequencing device connected to computer with SpotON Flow Cell inserted</td> <td><input type="checkbox"/> Desktop Agent set up (if applicable)</td> <td><input type="checkbox"/> Flow cell check complete</td> </tr> </table>		<input type="checkbox"/> Prepared library on ice	<input type="checkbox"/> Computer set up to run MinKNOW	<input type="checkbox"/> Hardware check complete	<input type="checkbox"/> Sequencing device connected to computer with SpotON Flow Cell inserted	<input type="checkbox"/> Desktop Agent set up (if applicable)	<input type="checkbox"/> Flow cell check complete
<input type="checkbox"/> Prepared library on ice	<input type="checkbox"/> Computer set up to run MinKNOW	<input type="checkbox"/> Hardware check complete					
<input type="checkbox"/> Sequencing device connected to computer with SpotON Flow Cell inserted	<input type="checkbox"/> Desktop Agent set up (if applicable)	<input type="checkbox"/> Flow cell check complete					

# 1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Priming and loading the SpotON flow cell</b></p>	
<p><input type="checkbox"/> Thaw one tube of Running Buffer with Fuel Mix (RBF) and one tube of Library Loading Beads (LBB) at RT and mix thoroughly by vortexing. Place on ice.</p> <p><input type="checkbox"/> Flip back the MinION lid and slide the priming port cover clockwise so that the priming port is visible.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Care must be taken when drawing back buffer from the flow cell. The array of pores must be covered by buffer at all times. Removing more than 20-30 <math>\mu</math>l risks damaging the pores in the array.</p>	
<p>After opening the priming port, check for small bubble under the cover. Draw back a small volume to remove any bubble (a few <math>\mu</math>l):</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 <math>\mu</math>l</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 <math>\mu</math>l, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p>Prepare the flow cell priming mix in a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 576 <math>\mu</math>l RBF</li> <li><input type="checkbox"/> 624 <math>\mu</math>l Nuclease-free water</li> </ul> <p><input type="checkbox"/> Load 800 <math>\mu</math>l of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes.</p> <p><input type="checkbox"/> Thoroughly mix the contents of the RBF and LLB tubes by pipetting.</p> <p>Prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 35.0 <math>\mu</math>l RBF</li> <li><input type="checkbox"/> 25.5 <math>\mu</math>l LLB</li> <li><input type="checkbox"/> 2.5 <math>\mu</math>l Nuclease-free water</li> <li><input type="checkbox"/> 12 <math>\mu</math>l DNA library</li> </ul> <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li> <li><input type="checkbox"/> Load 200 <math>\mu</math>l of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li> <li><input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading.</li> <li><input type="checkbox"/> Add 75 <math>\mu</math>l of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</li> <li><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION lid.</li> </ul>	
<p><b>Starting a sequencing run</b></p>	
<p><input type="checkbox"/> Double-click the MinKNOW icon located on the desktop to open the MinKNOW GUI.</p>	

Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> If your MinION was disconnected from the computer, plug it back in.</p> <p><input type="checkbox"/> Choose the flow cell type from the selector box. Then mark the flow cell as "Selected".</p> <p><input type="checkbox"/> Click the "New Experiment" button at the bottom left of the GUI.</p> <p>On the New experiment popup screen, select the running parameters for your experiment from the individual tabs.</p> <p><input type="checkbox"/> Output settings - FASTQ: The number of basecalls that MinKNOW will write in a single file. By default this is set to 4000</p> <p><input type="checkbox"/> Output settings - FAST5: The number of files that MinKNOW will write to a single folder. By default this is set to 4000</p> <p><input type="checkbox"/> Click "Start run".</p> <p>Allow the script to run to completion.</p> <p><input type="checkbox"/> The MinKNOW Experiment page will indicate the progression of the script; this can be accessed through the "Experiment" tab that will appear at the top right of the screen</p> <p><input type="checkbox"/> Monitor messages in the Message panel in the MinKNOW GUI</p> <p>The basecalled read files are stored in :data\reads</p>	
<p><b>Progression of MinKNOW protocol script</b></p>	
<p>The running experiment screen</p> <p>Experiment summary information</p> <p>Check the number of active pores reported in the MUX scan are similar (within 10-15%) to those reported at the end of the Flow Cell Check</p> <p><input type="checkbox"/> If there is a significant reduction in the numbers, restart MinKNOW.</p> <p><input type="checkbox"/> If the numbers are still significantly different, close down the host computer and reboot.</p> <p><input type="checkbox"/> When the numbers are similar to those reported at the end of the Flow Cell Check, restart the experiment on the Connection page. There is no need to load any additional library after restart.</p> <p><input type="checkbox"/> Stopping the experiment is achieved by clicking "Stop run" button at the top of the screen.</p> <p>Data acquisition will stop, but the software will continue basecalling unless the user clicks the "Stop basecalling" button.</p> <p><input type="checkbox"/> Check that the temperature has reached 34° C.</p> <p>Check pore occupancy in the channel panel at the top of the experimental view.</p> <p><input type="checkbox"/> A good library will be indicated by a higher proportion of light green channels in Sequencing than are in Pore. The combination of Sequencing and Pore indicates the number of active pores at any point in time. A low proportion of Sequencing channels will reduce the throughput of the run.</p> <p><input type="checkbox"/> Recovering indicates channels that may become available for sequencing again. A high proportion of this may indicate additional clean up steps are required during your library preparation.</p> <p><input type="checkbox"/> Inactive indicates channels that are no longer available for sequencing. A high proportion of these as soon as the run begins may indicate an osmotic imbalance.</p> <p><input type="checkbox"/> Unclassified are channels that have not yet been assigned one of the above classifications</p> <p><input type="checkbox"/> Monitor the pore occupancy</p>	



1D PCR barcoding genomic DNA (SQK-LSK108)

Version: PBGE\_9007\_v108\_revU\_18Oct2016  
 Last update: 03/09/2018



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Duty time plots</p> <p><input type="checkbox"/> Monitor the development of the read length histogram.</p> <p>Cumulative throughput</p> <p><input type="checkbox"/> the number of reads that have been sequenced and basecalled; and whether the reads have passed of failed the quality filters</p> <p>Trace viewer</p>	
<p><b>Onward analysis of MinKNOW basecalled data</b></p> <p><input type="checkbox"/> Open the Desktop Agent using the desktop shortcut.</p> <p><input type="checkbox"/> Click on the New Workflow tab in the Desktop Agent and select the FASTQ barcoding workflow.</p> <p>Select the workflow parameters.</p> <p><input type="checkbox"/> Select the quality score cut-off (this defaults to 7 unless changed)</p> <p><input type="checkbox"/> Select "Yes" in answer to "Detect barcode?"</p> <p><input type="checkbox"/> If you are working with human data, please tick "Yes" in answer to "Is the data you are about to upload a whole or partial human genome?", and confirm that you have consent from the subject to upload the data.</p> <p><input type="checkbox"/> Check the correct settings are selected in the Desktop Agent.</p> <p><input type="checkbox"/> Click "Start Run" to start data analysis.</p> <p><input type="checkbox"/> Follow the progression of upload and download of read files in the Desktop Agent.</p> <p>Click on VIEW REPORT.</p> <p><input type="checkbox"/> Click on VIEW REPORT to navigate to the Metrichor website, this can be done at any point during data exchange</p> <p><input type="checkbox"/> Return to the Desktop Agent to see progression of the exchange</p>	
<p><b>Close down MinKNOW and the Desktop Agent</b></p> <p><input type="checkbox"/> Quit Desktop Agent using the close x.</p> <p><input type="checkbox"/> Quit MinKNOW by closing down the web GUI.</p> <p><input type="checkbox"/> Disconnect the MinION.</p>	
<p><b>Prepare the flow cell for re-use or return to Oxford Nanopore.</b></p> <p><input type="checkbox"/> If you would like to reuse the flow cell, follow the Wash Kit instructions and store the washed flow cell at 2-8 °C, OR</p> <p><input type="checkbox"/> Follow the returns procedure by washing out the MinION Flow Cell ready to send back to Oxford Nanopore.</p>	