

Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 100 ng high molecular weight genomic DNA	<input type="checkbox"/> Agencourt AMPure XP beads	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
<input type="checkbox"/> PCR Barcoding Kit (SQK-PBK004)	<input type="checkbox"/> NEBNext End repair / dA-tailing Module (E7546)	<input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (M0367)	<input type="checkbox"/> Microfuge
	<input type="checkbox"/> Covaris g-TUBE	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Ice bucket with ice
	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Timer
	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937)	<input type="checkbox"/> Thermal cycler
	<input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> LongAmp Taq 2X Master Mix (e.g. NEB M0287)	
	<input type="checkbox"/> 10 mM Tris-HCl pH 8.0 with 50 mM NaCl	
	<input type="checkbox"/> (optional) Exonuclease I (NEB, M0293)	

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>DNA fragmentation</p> <p>OPTIONAL</p> <p>Prepare the DNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 100 ng genomic DNA into a DNA LoBind tube <input type="checkbox"/> Adjust the volume to 50 µl with Nuclease-free water <input type="checkbox"/> Mix thoroughly by inversion avoiding unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <p><input type="checkbox"/> Transfer 100 ng genomic DNA in 50 µl to the Covaris g-TUBE.</p> <p>Spin the g-TUBE for 1 minute at RT (Eppendorf 5424; 6000 rpm for 8 kb fragments).</p> <ul style="list-style-type: none"> <input type="checkbox"/> Spin the g-TUBE for 1 minute <input type="checkbox"/> Remove and check all the DNA has passed through the g-TUBE <input type="checkbox"/> If DNA remains in the upper chamber, spin again for 1 minute at the same speed 	

PCR Barcoding Kit (SQK-PBK004)

Version: PBK_9073_v1_revJ_14Aug2019
 Last update: 27/09/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Invert the g-TUBE and spin again for 1 minute to collect the fragmented DNA.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Remove g-TUBE, invert the tube and replace into the centrifuge <input type="checkbox"/> Spin the g-TUBE for 1 minute <input type="checkbox"/> Remove and check the DNA has passed into the lower chamber <input type="checkbox"/> If DNA remains in the upper chamber, spin again for 1 minute <input type="checkbox"/> Remove g-TUBE <p><input type="checkbox"/> Transfer the 50 µl fragmented DNA to a clean 1.5 ml Eppendorf DNA LoBind tube.</p>	
<p>100 ng fragmented DNA in 50 µl is taken into the next step.</p>	
<p>End-prep</p> <p>Perform end repair and dA-tailing of fragmented DNA as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 50 µl 100 ng fragmented DNA <input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer <input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Transfer the sample to a 0.2 ml PCR tube, and incubate for 5 minutes at 20° C and 5 minutes at 65° C using the thermal cycler.</p> <p><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</p> <p><input type="checkbox"/> Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual 70% ethanol. Briefly allow to dry.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 16 µl Nuclease-free water. Incubate for 2 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 16 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer.</p>	

PCR Barcoding Kit (SQK-PBK004)

Version: PBK_9073_v1_revJ_14Aug2019
 Last update: 27/09/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>PCR adapters ligation and amplification</p> <p>Add the reagents in the order given below, mixing by flicking the tube between each sequential addition:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 15 µl End-prepped DNA <input type="checkbox"/> 10 µl Barcode Adapters (BCA) <input type="checkbox"/> 25 µl Blunt/TA Ligase Master Mix <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. <input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing. <input type="checkbox"/> Add 20 µl of resuspended AMPure XP beads to the reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 21 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <p>Remove and retain 21 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Remove and retain the eluate which contains the DNA library in a clean 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Dispose of the pelleted beads <ul style="list-style-type: none"> <input type="checkbox"/> Quantify 1 µl of adapted DNA using a Qubit fluorometer. <input type="checkbox"/> Calculate how much DNA to take forward into the PCR step for a final DNA concentration of 0.2 ng/µl in a 50 µl reaction. <p>Set up the adapted DNA PCR as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Adapter ligated DNA, diluted x µl 0.2 ng/µl <input type="checkbox"/> Nuclease-free water 24-x µl <input type="checkbox"/> Barcode Primers (LWB 01-12, at 10 µM) 1 µl <input type="checkbox"/> LongAmp Taq 2x Master Mix 25 µl <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. 	

PCR Barcoding Kit (SQK-PBK004)

Version: PBK_9073_v1_revJ_14Aug2019
 Last update: 27/09/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Amplify using the following cycling conditions:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Initial denaturation 3 mins @ 95 °C (1 cycle) <input type="checkbox"/> Denaturation 15 secs @ 95 °C (14 (b) cycles) <input type="checkbox"/> Annealing 15 secs (a) @ 56 °C (a) (14 (b) cycles) <input type="checkbox"/> Extension 50 secs/kb @ 65 °C (c) (14 (b) cycles) <input type="checkbox"/> Final extension 6 mins @ 65 °C (1 cycle) <input type="checkbox"/> Hold @ 4 °C <ul style="list-style-type: none"> <input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing. <input type="checkbox"/> Add 30 µl of resuspended AMPure XP beads to the reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 10 µl of 10 mM Tris.HCl pH 8.0 with 50 mM NaCl. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <p>Remove and retain 10 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Remove and retain the eluate which contains the DNA in a clean 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Dispose of the pelleted beads <ul style="list-style-type: none"> <input type="checkbox"/> Quantify 1 µl of adapted DNA using a Qubit fluorometer. <input type="checkbox"/> Pool all barcoded libraries in the desired ratios to a total of 50-100 fmoles in 10 µl of 10 mM Tris-HCl pH 8.0 with 50 mM NaCl. 	
<p>Rapid adapter ligation</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Add 1 µl RAP to the 10 µl amplified DNA library. <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 5 minutes at RT. 	
<p>The prepared library is used for loading into the MinION flow cell. Store the library on ice until ready to load.</p>	

PCR Barcoding Kit (SQK-PBK004)

Version: PBK_9073_v1_revJ_14Aug2019
 Last update: 27/09/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Priming and loading the SpotON flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT before placing the tubes on ice as soon as thawing is complete. <input type="checkbox"/> Mix the Sequencing Buffer (SQB) and Flush Buffer (FB) tubes by vortexing, spin down and return to ice. <input type="checkbox"/> Spin down the Flush Tether (FLT) tube, mix by pipetting, and return to ice. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Please note that the Sequencing Tether (SQT) tube will NOT be used in this protocol. It is provided in the kit for potential future product compatibility. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Open the lid of the nanopore sequencing device and slide the flow cell's priming port cover clockwise so that the priming port is visible. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Care must be taken when drawing back buffer from the flow cell. The array of pores must be covered by buffer at all times. Removing more than 20-30 µl risks a loss of sequencing channels. 	
<p>After opening the priming port, check for small bubble under the cover. Draw back a small volume to remove any bubble (a few µl):</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip <ul style="list-style-type: none"> <input type="checkbox"/> Prepare the flow cell priming mix: add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by pipetting up and down. <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes. <input type="checkbox"/> Thoroughly mix the contents of the SQB and LB tubes by pipetting. <p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 34 µl Sequencing Buffer (SQB) <input type="checkbox"/> 25.5 µl Loading Beads (LB), mixed immediately before use <input type="checkbox"/> 4.5 µl Nuclease-free water <input type="checkbox"/> 11 µl DNA library 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> The Loading Beads (LB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use. 	
<p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. 	

PCR Barcoding Kit (SQK-PBK004)

Version: PBK_9073_v1_revJ_14Aug2019
Last update: 27/09/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Add 75 µl of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION lid.	
Ending the experiment	
<input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Wash Kit instructions and store the washed flow cell at 2-8 °C, OR <input type="checkbox"/> Follow the returns procedure by washing out the flow cell ready to send back to Oxford Nanopore.	
IMPORTANT <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.	