

PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Version: PCB_9155_v111_revL_18May2022
 Last update: 07/03/2024



Flow Cell Number:

DNA Samples:

Before start checklist

Materials

- 4 ng enriched RNA (Poly(A)+ RNA or ribodepleted) or 200 ng total RNA
- PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Consumables

- Agencourt RNAClean XP beads (Beckman Coulter™, cat # A63987)
- Agencourt AMPure XP beads (Beckman Coulter™ cat # A63881)
- Lambda Exonuclease (NEB, Cat # M0262L)
- NEBNext® Quick Ligation Reaction Buffer (NEB, B6058)
- T4 DNA Ligase 2M U/ml (NEB, cat # M0202T/M)
- 10 mM dNTP solution (e.g. NEB N0447)
- LongAmp Hot Start Taq 2X Master Mix (NEB, M0533S)
- Maxima H Minus Reverse Transcriptase (200 U/μl) with 5x RT Buffer (ThermoFisher, cat # EP0751)
- RNaseOUT™, 40 U/μl (Life Technologies, cat # 10777019)
- USER (Uracil-Specific Excision Reagent) Enzyme (NEB, cat # M5505L)
- Exonuclease I (NEB, Cat # M0293)
- Nuclease-free water (e.g. ThermoFisher, AM9937)
- Freshly prepared 70% ethanol in nuclease-free water
- 1.5 ml Eppendorf DNA LoBind tubes
- 0.2 ml thin-walled PCR tubes
- Qubit RNA HS Assay Kit (ThermoFisher, cat # Q32852)
- Qubit dsDNA HS Assay Kit (ThermoFisher, cat # Q32851)
- Qubit™ Assay Tubes (Invitrogen, Q32856)

Equipment

- Hula mixer (gentle rotator mixer)
- Magnetic rack, suitable for 1.5 ml Eppendorf tubes
- Microfuge
- Vortex mixer
- Thermal cycler
- Ice bucket with ice
- Timer
- Qubit fluorometer (or equivalent for QC check)
- Agilent Bioanalyzer (or equivalent)
- Pipettes and pipette tips P2, P10, P20, P100, P200, P1000

INSTRUCTIONS

NOTES/OBSERVATIONS

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Reverse transcription and strand-switching</p> <p>Thaw the following reagents, then spin down briefly using a microfuge and mix as indicated in the table below. Then place the reagents on ice.</p> <ul style="list-style-type: none"> <input type="checkbox"/> cDNA RT Adapter (CRTA): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> Annealing Buffer (AB): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> Short Fragment Buffer (SFB): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> RT Primer (RTP): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> Strand Switching Primer II (SSPII): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> NEBNext® Quick Ligation Reaction Buffer: thaw at RT, briefly spin down, mix by vortexing <input type="checkbox"/> T4 DNA Ligase 2M U/ml: not frozen, briefly spin down, mix well by pipetting <input type="checkbox"/> RNaseOUT: not frozen, briefly spin down, mix well by pipetting <input type="checkbox"/> Lambda Exonuclease: not frozen, briefly spin down, mix well by pipetting <input type="checkbox"/> Uracil-Specific Excision Reagent (USER): not frozen, briefly spin down, mix well by pipetting <input type="checkbox"/> 10 mM dNTP solution: thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> Maxima H Minus Reverse Transcriptase: not frozen, briefly spin down, mix well by pipetting <input type="checkbox"/> Maxima H Minus 5x RT Buffer: thaw at RT, briefly spin down, mix by vortexing 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> It is important that the NEBNext Quick Ligation Reaction Buffer is mixed well by vortexing. 	
<p>For each sample, prepare the RNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 4 ng Poly(A)+ RNA, or 200 ng total RNA into a 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Adjust the volume up to 10 µl with Nuclease-free water <input type="checkbox"/> Mix by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <p>Prepare the following in a 0.2 ml PCR tube per sample:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 10 µl RNA <input type="checkbox"/> 1 µl cDNA RT Adapter (CRTA) <input type="checkbox"/> 1 µl Annealing Buffer (AB) <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tubes, and spin down. <input type="checkbox"/> Incubate the reactions in the thermal cycler at 60°C for 5 mins, then cool for 10 minutes at RT. <p>To each of the 0.2 ml PCR tubes containing you RNA sample(s), add the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 12 µl RNA sample (from previous step) <input type="checkbox"/> 3.6 µl NEBNext® Quick Ligation Reaction Buffer <input type="checkbox"/> 1.4 µl T4 DNA Ligase 2M U/ml <input type="checkbox"/> 1 µl RNaseOUT <ul style="list-style-type: none"> <input type="checkbox"/> Ensure the components are thoroughly mixed by flicking the tubes and spin down. <input type="checkbox"/> Incubate for 10 minutes at RT. 	

PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Version: PCB_9155_v111_revL_18May2022
 Last update: 07/03/2024



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>To each of the 0.2 ml PCR tubes, add the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 18 µl RNA sample (from previous step) <input type="checkbox"/> 1 µl Lambda Exonuclease <input type="checkbox"/> 1 µl USER (Uracil-Specific Excision Reagent) <p><input type="checkbox"/> Ensure the components are thoroughly mixed by flicking the tubes and spin down.</p> <p><input type="checkbox"/> Incubate for 15 minutes at 37°C in the thermal cycler.</p> <p><input type="checkbox"/> Transfer each sample to clean 1.5 ml Eppendorf DNA LoBind tubes.</p> <p><input type="checkbox"/> Resuspend the RNase-free XP beads by vortexing.</p> <p><input type="checkbox"/> Add 36 µl of resuspended RNase-free XP beads to each reaction and mix gently by flicking the tubes.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Spin down the samples and pellet on a magnet. Keep the tubes on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tubes on the magnet and wash the beads with 100 µl of Short Fragment Buffer (SFB) without disturbing the pellet. Remove the SFB using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tubes back on the magnet. Pipette off any residual buffer. Briefly allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tubes from the magnetic rack and resuspend each pellet in 12 µl of Nuclease-free water.</p> <p><input type="checkbox"/> Incubate at RT for 10 minutes.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 12 µl of eluate into a clean 0.2 ml thin-walled PCR tube per sample.</p> <p>To each of the 0.2 ml PCR tubes, add the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 12 µl Eluted sample (from previous step) <input type="checkbox"/> 1 µl RT Primer (RTP) <input type="checkbox"/> 1 µl dNTPs (10 mM) <p><input type="checkbox"/> Ensure the components are thoroughly mixed by flicking the tubes and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 15 minutes at RT.</p> <p>To each of the 0.2 ml PCR tubes, add the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 14 µl RT primed sample (from previous step) <input type="checkbox"/> 4.5 µl Maxima H Minus 5x RT Buffer <input type="checkbox"/> 1 µl RNaseOUT <input type="checkbox"/> 2 µl Strand Switching Primer II (SSP II) 	

PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Version: PCB_9155_v111_revL_18May2022
 Last update: 07/03/2024



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Mix gently by flicking the tubes, and spin down.</p> <p><input type="checkbox"/> Incubate at 42°C for 2 minutes in the thermal cycler.</p> <p><input type="checkbox"/> Add 1 µl of Maxima H Minus Reverse Transcriptase to each tube. The total volume will be 22.5 µl per tube.</p> <p><input type="checkbox"/> Mix gently by flicking the tubes, and spin down.</p> <p>Incubate using the following protocol using a thermal cycler:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Reverse transcription and strand-switching 90 mins @ 42°C (1 cycle) <input type="checkbox"/> Heat inactivation 5 mins @ 85°C (1 cycle) <input type="checkbox"/> Hold @ 4°C 	
<p>Take your samples forward into the next step. However, at this point it is also possible to store the sample at -20°C overnight.</p>	
<p>Selecting for full-length transcripts by PCR</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> This kit enables multiplexing of up to 24 samples. The default method allows you to perform one 25 µl PCR reaction per sample. If multiplexing two or three samples, however, two separate PCR reactions per sample should be performed; if running just one sample, four separate PCR reactions should be performed as per the cDNA-PCR Sequencing Kit protocol (SQK-PCS111). These recommendations aim to ensure that enough PCR product is generated for optimal flow cell performance.</p>	
<p>Thaw the following reagents, then spin down briefly using a microfuge and mix as indicated in the table below. Then place the reagents on ice.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Barcode Primers (BP01 - BP24): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> Elution Buffer (EB): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> LongAmp Hot Start Taq 2X Master Mix: thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> Exonuclease I: not frozen, briefly spin down, mix well by pipetting <p><input type="checkbox"/> Spin down the reverse-transcribed RNA samples.</p> <p><input type="checkbox"/> Prepare a separate 0.2 ml PCR tube for each sample and add 5 µl of reverse-transcribed RNA per tube.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> Only 5 µl of the reverse-transcribed sample is to be taken forward. Do NOT use all the 22.5 µl of the reverse transcription reaction in a single PCR reaction.</p>	
<p>In each of the 0.2 ml PCR tubes containing reverse-transcribed RNA sample, prepare the following reaction at RT:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 5 µl Reverse-transcribed sample (from previous step) <input type="checkbox"/> 0.75 µl Unique Barcode Primer (BP01-24) <input type="checkbox"/> 6.75 µl Nuclease-free water <input type="checkbox"/> 12.5 µl 2x LongAmp Hot Start Taq Master Mix <p><input type="checkbox"/> Mix gently by pipetting.</p>	

PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Version: PCB_9155_v111_revL_18May2022
 Last update: 07/03/2024



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Amplify using the following cycling conditions.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Initial denaturation 30 secs @ 95°C (1 cycle) <input type="checkbox"/> Denaturation 15 secs @ 95°C (10-18* cycles) <input type="checkbox"/> Annealing 15 secs @ 62°C (10-18* cycles) <input type="checkbox"/> Extension 60 secs per kb @ 65°C (10-18* cycles) <input type="checkbox"/> Final extension 6 mins @ 65°C (1 cycle) <input type="checkbox"/> Hold @ 4°C <p><input type="checkbox"/> Add 1 µl Exonuclease I directly to each PCR tube. Mix by pipetting.</p> <p><input type="checkbox"/> Incubate the reactions at 37°C for 15 minutes, followed by 80°C for 15 minutes in the thermal cycler.</p> <p><input type="checkbox"/> Transfer each sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Add 20 µl of resuspended AMPure XP beads to each 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 5 ml of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the samples and pellet on a magnet. Keep the tubes on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tubes on the magnet and wash the beads with 200 µl of freshly-prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tubes back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellets to the point of cracking.</p> <p><input type="checkbox"/> Remove the tubes from the magnetic rack and resuspend each pellet in 12 µl of Elution Buffer (EB).</p> <p><input type="checkbox"/> Incubate at RT for 10 minutes.</p> <p><input type="checkbox"/> Pellet the beads on the magnet until the eluate is clear and colourless.</p> <p>Remove and retain 12 µl of each eluate into a separate clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Remove and retain the eluate which contains the cDNA library in a clean 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Dispose of the pelleted beads <p><input type="checkbox"/> For each sample, analyse 1 µl of the amplified cDNA for size, quantity and quality using a Qubit fluorometer and Agilent Bioanalyzer (or equivalent) for a QC check.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> Sometimes a high-molecular weight product is visible in the wells of the gel when the PCR products are run, instead of the expected smear. These libraries are typically associated with poor sequencing performance. We have found that repeating the PCR with fewer cycles can remedy this.</p>	

PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Version: PCB_9155_v111_revL_18May2022
 Last update: 07/03/2024



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Pool together equimolar samples of the amplified cDNA barcoded samples to a total of 15 – 25 fmols and make the volume up to 23 µl in Elution Buffer (EB).</p> <p><input type="checkbox"/> Please check the Mass to Molarity table in the protocol</p>	
<p>Adapter addition</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> The Rapid Adapter T (RAP T) used in this kit and protocol is not interchangeable with other sequencing adapters.</p> <p><input type="checkbox"/> Spin down the Rapid Adapter T (RAP T) and place on ice.</p> <p><input type="checkbox"/> Thaw the RAP Dilution Buffer (RDB) or Adapter Buffer (ADB) at RT, spin down briefly using a microfuge and mix by pipetting before storing on ice.</p> <p>In a fresh 1.5 ml Eppendorf DNA LoBind tube, dilute Rapid Adapter T (RAP T):</p> <p><input type="checkbox"/> 1.2 µl Rapid Adapter T (RAP T)</p> <p><input type="checkbox"/> 6.8 µl RAP Dilution Buffer (RDB) or Adapter Buffer (ADB)</p> <p><input type="checkbox"/> Mix well by pipetting and spin down.</p> <p><input type="checkbox"/> Add 1 µl of the diluted Rapid Adapter T (RAP T) to the amplified cDNA library, making the final volume up 24 µl.</p> <p><input type="checkbox"/> Mix well by pipetting and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 5 minutes at RT.</p>	
<p>The prepared library is used for loading onto the flow cell. Store the library on ice until ready to load.</p>	
<p>Priming and loading the flow cell</p>	
<p>Using the Loading Solution</p> <p><input type="checkbox"/> Thaw the Sequencing Buffer II (SBII), Loading Beads II (LBII) or Loading Solution (LS, if using), Flush Tether (FLT) and Flush Buffer (FB) at RT before mixing the reagents by vortexing, and spin down the SBII and FLT at RT.</p> <p><input type="checkbox"/> To prepare the flow cell priming mix, add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> After taking flow cells out of the fridge, wait 20 minutes before inserting the flow cell into the PromethION for the flow cell to come to RT. Condensation can form on the flow cell in humid environments. Inspect the gold connector pins on the top and underside of the flow cell for condensation and wipe off with a lint-free wipe if any is observed. Ensure the heat pad (black pad) is present on the underside of the flow cell.</p>	
<p>For PromethION 2 Solo, load the flow cell(s) as follows:</p> <p><input type="checkbox"/> Place the flow cell flat on the metal plate.</p> <p><input type="checkbox"/> Slide the flow cell into the docking port until the gold pins or green board cannot be seen.</p>	

PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Version: PCB_9155_v111_revL_18May2022
 Last update: 07/03/2024



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>For the PromethION 24/48, load the flow cell(s) into the docking ports:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Line up the flow cell with the connector horizontally and vertically before smoothly inserting into position. <input type="checkbox"/> Press down firmly onto the flow cell and ensure the latch engages and clicks into place. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Insertion of the flow cells at the wrong angle can cause damage to the pins on the PromethION and affect your sequencing results. If you find the pins on a PromethION position are damaged, please contact support@nanoporetech.com for assistance. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Slide the inlet port cover clockwise to open. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores. 	
<p>After opening the inlet port, draw back a small volume to remove any air bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette tip to 200 µl. <input type="checkbox"/> Insert the tip into the inlet port. <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you see a small volume of buffer entering the pipette tip. <ul style="list-style-type: none"> <input type="checkbox"/> Load 500 µl of the priming mix into the flow cell via the inlet port, avoiding the introduction of air bubbles. Wait five minutes. During this time, prepare the library for loading using the next steps in the protocol. <input type="checkbox"/> Thoroughly mix the contents of the Loading Beads II (LBII) by pipetting. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> The Loading Beads II (LBII) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use. 	
<p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 75 µl Sequencing Buffer II (SBII) <input type="checkbox"/> 51 µl Loading Beads II (LBII) thoroughly mixed before use, or Loading Solution (LS), if using <input type="checkbox"/> 24 µl DNA library <ul style="list-style-type: none"> <input type="checkbox"/> Complete the flow cell priming by slowly loading 500 µl of the priming mix into the inlet port. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Using a P1000, insert the pipette tip into the inlet port and add 150 µl of library. <input type="checkbox"/> Close the valve to seal the inlet port. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output. 	
<p>If the light shield has been removed from the flow cell, install the light shield as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Align the inlet port cut out of the light shield with the inlet port cover on the flow cell. The leading edge of the light shield should sit above the flow cell ID. <input type="checkbox"/> Firmly press the light shield around the inlet port cover. The inlet port clip will click into place underneath the inlet port cover. 	

PCR-cDNA Barcoding Kit (SQK-PCB111.24)

Version: PCB_9155_v111_revL_18May2022
Last update: 07/03/2024



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
Close the PromethION lid when ready to start a sequencing run on MinKNOW.	
Flow cell reuse and returns	
<input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C. <input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore.	
IMPORTANT <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.	